ACE (Angiotensin converting Enzyme)

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5th Ed. AACC Press (2000).

3. CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.

4. Maguire G.A. et al., Ann. Clin. Biochem. 22, 204 (1985).

5. SIGMA Diagnostics IFU 305-UV, rel. 10/1994

BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test,

ACE (Angiotensin converting Enzyme)

Angiotensin converting Enzyme (ACE) QUANTITATIVE KINETIC COLORIMETRIC ASSAY
FAPGG (N-[3-(2-furyl)acryloyl]-L-phenylalanyl-glycylglycine) METHOD on SERUM and PLASMA with MINDRAY BS-300 Analyzer

500 tests (10 x 10 mL)

NACE8865

INTENDED USE

Angiotensin-converting enzyme (ACE) is an exopeptidase. It is a circulating enzyme that participates in the body's renin-angiotensin-aldosterone system (RAAS); this one is a hormone system that regulates blood pressure and water fluid balance.

When blood volume is low, JG cells or granular cells in the kidneys secrete renin directly into circulation; plasma renin then converts angiotensinogen released by the liver into angiotensin I. Angiotensin I is subsequently converted to angiotensin II by ACE (angiotensin converting enzyme), enzyme found in the lungs.

Angiotensin II is a potent vaso-active peptide that causes blood vessels to constrict, giving increased blood pressure. Angiotensin II also stimulates the secretion of the aldosterone, a hormone from adrenal cortex.

This one causes the tubules of the kidneys to increase the reabsorption of sodium and water into the blood; this increases the volume of fluid in the body, which also increases blood pressure.

If RAAS is too active, blood pressure will be too high; many drugs are able to interrupt different steps in this system to lower blood pressure.

ACE degrades bradykinin, a potent vasodilator, and other vasoactive peptides.

Inhibition of ACE by ACE-inhibitors results in the decreased formation of angiotensin II and decreased catabolism of bradykinin, leading to systematic dilation of the arteries and veins and a decrease in blood pressure.

In addition, inhibiting angiotensin II formation diminishes angiotensin II-mediated aldosterone secretion, leading to a decrease in water and sodium reabsorption and a reduction in extracellular volume.

An increase of ACE values happens for sarcoidosis, leprosy, arterial hypertension, rheumatoid arthritis, acute and chronic bronchitis, **hyperthyroidism**, cirrhosis of the liver.

PRINCIPLE

The ACE present in the serum catalyzes the hydrolysis of the FAPGG, forming furylacryloylphenylalanine (FAP).

The decrease of absorbance in the unit time at 340 nm is proportional at the activity of the ACE in the sample.

Using our Calibrator **REF ACECAL8866**, it is possible to prepare a Calibration Curve to refer.

Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample.

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change

the result.

3. DO NOT mix Reagents from different Production lots.
4. For concentration of ACE higher than 150 U/L, dilute the sample 1:4 with saline solution, repeat the determination and multiply the result by 4.

5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of ACE (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

G) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

REAGENTS

Components of the kit:

R1 - BUFFER

Good Buffer >20 mmol/L

R2 - FAPGG

FAPGG > 0.25 mmol/LL

NaN3 < 0.1%

REF NACE8865

2 x 50 mL

10 x 10 mL

STABILITY: the Reagents are stable up to the expiry date

Mascia Brunelli S.p.a.

mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- ACE CALIBRATOR Iyo
- ACE NORMAL CONTROL Iyo
- ACE ELEVATED CONTROL Iyo

REF	ACECAL8866
REF	ANOR8867
REF	AEL8868

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - FAPGG** with 10 mL of **R1 - BUFFER**. Mix gently and let the reagent reach the room temperature before use. Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination. An incompetent handling relieves us from any responsibility.

STABILITY OF THE WORKING REAGENT

The WORKING REAGENT is stable 4 weeks at 2-8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

- Normal laboratory equipment.
- Micropipettes to deliver from 3 to 1000 µL.
- Disposable micropipettes tips .
- Transparent glass tubes for sample dilution.
- Saline solution, distilled water, Calibrators and Controls.
- Spectrophotometer or automatic analyzer for Clinical Chemistry.

SAMPLES

- Not haemolysed fresh serum.
- Not haemolysed fresh plasma plasma with heparin. EDTA cannot be used!

Samples collection in compliance with CLSI (NCCLS) (see References 3).

The sample can be stored at 2-8°C, up to 6 days. (see References 1).

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

ANALYTICAL SPECTROPHOTOMETER PROCEDURE ON MANUAL

- Wavelength: 340 nm
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: kinetic
- Reaction: 5 + 5 minutes
- Sample/reagent: 1/10

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, ST: Calibrator, S: Sample

	S	ST
--	---	----

WORKING REAGENT	1000 µL	1000 µL
Sample	100 µL	----
Calibrator	----	100 µL

Mix wells, incubate for 5 minutes at 37°C.

Read the absorbance of standard (Ast1) and sample (As1). Exactly after 5 minutes at 37°C read again standard (Ast2) and sample (As2).

Determine the diff. of absorbance for sample and calibrator:

$$\Delta A_s = A_{1s} - A_{2s}$$

$$\Delta A_{st} = A_{1st} - A_{2st}$$

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analysers could be completely different by what has been developed as manual determination.

REFERENCE VALUES (see References 5)

Normal Values ACE: 8 - 52 U/L.

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

CALCULATION

Insert the means found in the following formula:

$$(\Delta A_s / \Delta A_{st}) \times \text{Calibrator conc.} = \text{U/L of ACE.}$$

ANALYTICAL PERFORMANCES

(validate on MINDRAY BS300)

The performances of the Reagent **ACE (Angiotensin converting Enzyme)** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 150 U/L.

However, for ACE concentrations higher than 150 U/L, it is recommended to dilute the sample 1:4 with saline solution, test again and multiply the result x 4.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 2.6 U/L.

Interferences: see References 2.

Interference test criterion: recovery \pm 15% of initial value. No interference found on samples with:

- total bilirubin up to 20 mg/dL;
- haemoglobin up to 300 mg/dL,
- lipemia [Intralipid®] up to 1000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	39.9 \pm 1.2	1.5
Human 2	100.3 \pm 1.5	0.7

Mascia Brunelli S.p.a.

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	39.4 \pm 2.1	2.6
Human 2	99.4 \pm 2.8	1.4

Accuracy: a group of 20 sera has been tested using this procedure and using the reagent on different analyzers.

The comparison gave these results:

Linear regression equation $y = 0.9985x + 0.475$

Correlation coefficient $r = 0.9995$ $n = 20$

REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.4. Maguire G.A. et al., Ann. Clin. Biochem. 22, 204 (1985).
5. SIGMA Diagnostics IFU 305-UV, rel. 10/1994.