PYRUVATE

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Regressione lineare $y = 1.1252x - 0.034$

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BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.

PYRUVATE

PYRUVATE (PYRUVIC ACID) QUANTITATIVE UV ASSAY on WHOLE BLOOD
with MANUAL SPECTROPHOTOMETER HP8452A

100 tests (10 x 10 mL)

REF NAPY8825

INTENDED USE

Pyruvate (pyruvic acid) is a key intermediate in cellular metabolic pathways. Pyruvate is primarily derived from glucose via glycolysis. In the case that sufficient oxygen is available, pyruvate is converted by pyruvate dehydrogenase into acetyl-CoA, which enters the citric acid cycle (Krebs cycle), where it is metabolized to produce ATP.

Pyruvate can also be converted into carbohydrates via gluconeogenesis or into fatty acids through acetyl-CoA.

In the last years, pyruvate has been used as a dietary supplement to aid in weight loss. Pyruvate can improve exercise endurance capacity, effectively reduces cholesterol and serves as a potent antioxidant. Therefore the measurement of pyruvate concentration can give valuable information to the progress of specific biochemical reactions.

The lactate to pyruvate ratio reflects the redox state of the cell and describe the balance between NAD and NADH, which is dependent on the interconversion of lactate and pyruvate by lactate dehydrogenase (LDH).

This ratio has been used to assess the severity of circulatory failure as well as other disease states.

Abnormal increase of blood pyruvate levels are reported in a number of disorders, as shock, liver disease, congestive heart failure, neoplastic disorders, muscular dystrophy, diabetes mellitus, thiamine deficiency and metabolic disorders.

PRINCIPLE

The Lactate Dehydrogenase (LDH) catalyzes the reduction of pyruvate to lactate with simultaneous oxidation of NADH to NAD. The decrease of absorbance of NADH, is directly proportional to the pyruvate concentration in the sample.

Using the standard contained in the kit it is possible to prepare a Calibration Curve to refer.

Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample.

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of Pyruvate higher than 4.1 mg/dL (0.460 mmol/L), dilute the sample 1:2 with saline solution, repeat the determination and multiply the result by 2.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous

membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

A) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Pyruvate (see References 2).

B) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

C) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

D) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

E) The R1 - BUFFER is supplied in surplus.

F) The preparation of working reagent and the procedure for calibration are different from those for samples and controls.

G) Controls and samples will be deproteinized before the test.

H) R4 - CAL is a standard ready to use. Do not deproteinize!

I) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

J) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

REAGENTS

Components of the kit:

R1 - BUFFER	Good buffer > 400 mmol/L	REF NAPY8825 3 x 70 mL
R2 - NADH	LDH >250 KU/L	10 x 10 mL
R3 - LDH (READY TO USE)	NADH >0,1 mmol/L Activators, Stabilizers	1 x 5 mL
R4 - CAL (READY TO USE)	Solution of Pyruvate = 4.0 mg/dL (= 0.450 mmol/L) NaN3 < 0.1%	1 x 10 mL

STABILITY: the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- SUBSTRATE ELEVATED CONTROL Iyo
 - SUBSTRATE LOW CONTROL Iyo

REF OG3005
 REF OG3006

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

PREPARATION OF THE REAGENT FOR CALIBRATION (ONLY FOR CALIBRATION)

Dissolve a vial of **R2 - NADH** with 10 mL of **R1 - BUFFER**, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

STABILITY OF THE REAGENT FOR CALIBRATION (ONLY FOR CALIBRATION)

The **REAGENT FOR CALIBRATION** is stable 7 days at 2-8°C.

PREPARATION OF THE WORKING REAGENT (ONLY FOR CONTROLS AND SAMPLES)

Dissolve a vial of **R2 - NADH** with 5 mL of **R1 - BUFFER**, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

STABILITY OF THE WORKING REAGENT (ONLY FOR CONTROLS AND SAMPLES)

The **WORKING REAGENT** is stable 7 days at 2-8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Saline solution, distilled water, Controls.

Perchloric acid 70%.

Spectrophotometer.

SAMPLES

• Whole blood.

Samples collection in compliance with CLSI (NCCLS)

(see References 3).

For samples stability, see the following paragraph.

SAMPLES PRETREATMENT

Prepare protein-free supernatants as follows:

- draw **BLOOD** following standardized procedure: It is preferable the patient be in a fasting and resting state;

- into a centrifuge tube put 4.0 mL of cold 8% (w/v) perchloric acid (PCA) (prepared by diluting 7 mL of 70% (w/w) perchloric acid to 100 mL with distilled water);

- quick pipet 2.0 mL **BLOOD** into that tube and vortex mix for 30 seconds ;

- keep cold the blood-precipitate mixture again for 5 minutes, to assure complete protein precipitation;

- centrifuge for about 10 minutes x 1500 g ;

- the clear PCA supernatant is ready for use.

NOTE: - a second centrifugation of the supernatant may be necessary to obtain a clear protein-free solution;

- all the operation must require a short time; pyruvate in PCA is stable for a few weeks refrigerated.

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

ANALYTICAL PROCEDURE FOR CALIBRATION ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 5 + 5 minutes

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled:

Cuvet n°	Reagent for Calibration	Distilled water	R4 - CAL	BLOOD PYRUVATE (mg/dL (mmol/L))	
1	1000 µl	1900 µl	100 µl	0.600	0.068
2	1000 µl	1800 µl	200 µl	1.200	0.136
3	1000 µl	1700 µl	300 µl	1.800	0.204
4	1000 µl	1600 µl	400 µl	2.400	0.272
5	1000 µl	1500 µl	500 µl	3.000	0.340

Mix kindly and incubate for 5 minutes at 37°C.

Read for the first time the absorbance of sample (As1) against Water (INITIAL Abs). Immediately after this one, add to the 5 cuvetts:

R3 - LDH	50 µl
-----------------	--------------

Mix kindly. Exactly 5 minutes after the first reading, make the second reading of sample (As2) (FINAL Abs).

ANALYTICAL PROCEDURE FOR SAMPLES AND CONTROLS ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 5 + 5 minutes
- Sample/reagent: 4/1/1/0.1

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvet so labelled :

S: Sample:

	S
PCA supernatant fluid (from sample deproteinization)	2000 µl
R1 - BUFFER	500 µl

Mix kindly but surely: the mixture must be at the right pH before adding **WORKING REAGENT**. Incubate for 5 minutes x 37°C and add:

WORKING REAGENT	500 µl
------------------------	--------

Read for the first time the absorbance of sample (As1) against Water (INITIAL Abs). Immediately after this one, add:

R3 - LDH	50 µl
-----------------	-------

Mix kindly. Exactly 5 min. after the first reading, make the second

reading of sample (As2) (FINAL Abs).

REFERENCE VALUES (see References 1)

Normal Values Pyruvic Acid 0.3 - 0.9 mg/dL (0.034 - 0.102 mmol/L).

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

CALCULATION

Plot each value found on the Calibration Curve. The Calibration Curve has to be always repeated for each new lot of Reagent;

or

$(As1 - As2) \times 6.396 = \text{mg/dL of Pyruvic Acid}$

$(As1 - As2) \times 0.726 = \text{mmol/L of Pyruvic Acid}$

Pyruvic Acid in mmol/L = $0.1135 \times \text{mg/dL of Pyruvic Acid}$

ANALYTICAL PERFORMANCES

(validate on SPECTROPHOTOMETER HP8452A)

The performances of the Reagent **PYRUVATE** have been tested with a spectrophotometer HP8452A. The data, while representing the characteristics of the product, could be different for each laboratory and for different analysers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 4.1 mg/dL

(0.460 mmol/L). For concentrations ≥ 4.1 mg/dL (0.460 mmol/L) it is recommended to dilute the sample 1:2 with saline solution, test again and multiply the result x 2.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 0 mg/dL (0.00 mmol/L).

Interferences: see References 2.

Interference test criterion: recovery $\pm 10\%$ of initial value. No interference found on samples with:

- total bilirubin up to 10 mg/dL;
- haemoglobin up to 600 mg/dL,
- lipemia [Intralipid®] up to 4000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (mmol/L) $\pm 2s$	CV %
Human 1	0.18 \pm 0.02	5.9
Human 2	0.34 \pm 0.04	6.6

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (mmol/L) $\pm 2s$	CV %
Human 1	0.18 \pm 0.02	7.1
Human 2	0.34 \pm 0.02	2.4

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market.

The comparison gave these results:

Linear regression equation $y = 1.1252x - 0.034$

Correlation coefficient $r = 0.9908$ $n = 20$

REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).

2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

3. CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.