## PYRUVATE KINASE

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Codice Ramo **CND W01010199**

## PYRUVATE KINASE

QUANTITATIVE KINETIC UV DETERMINATION OF PYRUVATE KINASE (PK)  
in SERUM and ERYTHROCYTES on MINDRAY BS-300 Analyzer

265 tests (20 x 3 mL)

NCPK8831

### INTENDED USE

Pyruvate kinase (PK) catalyzes the second step of formation of energy in the glycolysis; a deficiency involves strict consequences here. The only scope of this enzyme is to catalyze the transfer of a phosphate group from phosphoenolpyruvate (PEP) to ADP, yielding one molecule of pyruvate and one molecule of ATP.

The activation of this enzyme needs magnesium and potassium. PK activity is regulated by:

-its own substrate, PEP and fructose 1,6-biphosphate, enhance enzymatic activity. Thus glycolysis is driven to operate faster with more substrate.

-ATP and alanine are negative inhibitor or modulator.

-acetyl-CoA seems to have negative inhibition; the potential inhibitory role of citrate is not well known.

Liver PK is also regulated indirectly by epinephrine (inhibition) and glucagon (activation), while muscle PK is not inhibited by epinephrine.

Genetic defects of this enzyme cause pyruvate kinase deficiency. In this condition, a lack of pyruvate kinase slows down the process of glycolysis. This effect is especially devastating in cells that lack mitochondria, because these cells must use anaerobic glycolysis as their sole source of energy because the Krebs cycle is not available. If an human being has red blood cells in a state of pyruvate kinase deficiency, rapidly they become deficient in ATP and can undergo to hemolytic anemia.

Sensitive variations of the PK are found in the tumors, in particular liver tumors.

### PRINCIPLE

PK in the sample catalyzes the reaction between ADP and PEP (phosphoenolpyruvate), forming Pyruvate.

The Lactate dehydrogenase (LDH) in the system catalyzes the reaction between Pyruvate and NADH, giving Lactic ac. and NAD. The decrease of absorbance in 10 min at 340 nm, due to the oxidation of NADH is proportional at the activity of the PK in the sample.

Using the calibrator contained in the kit it is possible to prepare a Calibration Curve to refer.

Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample.

### PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.

4. For concentration of Pyruvate Kinase higher than 90 U/L for serum and 900 mU/10<sup>9</sup> erythrocytes /mL, dilute the sample 1:10 with saline solution, repeat the determination and multiply the result by 10.

5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product.

Material safety data sheet should be supplied on request.

### ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Pyruvate Kinase (see References 2):

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) Whole blood sample will be pretreated before the test.

F) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

G) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

H) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

### REAGENTS

Components of the kit:

#### R1 - BUFFER

Good buffer >20 mmol/L

#### R2 - NADH, LDH

NADH >0.1 mmol/L  
PEP >0.1 mmol/L

REF NCPK8831

1 x 60 mL

20 x 3 mL

LDH	> 1000 U/L	
<b>R3 - ADP</b>		<b>3 x 2 mL</b>
ADP	>0.1 mmol/L	
<b>R4 - CAL</b>		<b>3 x 1 mL</b>
Human stabilized serum (Values of Pyruvate Kinase, see on the bottle)		
<b>R5 - DIL</b>		<b>1 x 4 mL</b>
Calibrator diluent NaN3 < 0.1%		

**STABILITY:** the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

#### AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- PIRUVATE KINASE (PK) CONTROL REF OGPCKO32

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

#### PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - NADH, LDH** with 3 mL of **R1 - BUFFER**, mix gently until dissolution.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

#### STABILITY OF THE WORKING REAGENT

The WORKING REAGENT is stable 3 days at 2-8°C.

#### PREPARATION OF THE WORKING STARTER

Dissolve a vial of **R3 - ADP** with 2 mL of DISTILLED WATER, mix gently until dissolution. Let the reagent reach the room temperature before use. Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

#### STABILITY OF THE WORKING STARTER

The WORKING STARTER is stable 2 weeks at 2-8°C.

#### PREPARATION OF CALIBRATOR

Dissolve a vial of **R4 - CAL** with 1 mL of **R5 - DIL**, mix gently until dissolution. Avoid foam. Let the reagent reach the room temperature before use. Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

#### STABILITY OF CALIBRATOR

The CALIBRATOR is stable:

- 4 hours at 2-8°C in the dark;

- 3 days at -20°C in the dark.

#### MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Saline solution, distilled water, Controls.

Spectrophotometer or automatic analyzer for Clinical Chemistry.

#### SAMPLES

• Not haemolysed and not lipemic fresh serum, without other treatments.

To test immediately. A 50% of activity may be lost after 24 h at 2-8°C.

• Whole blood collected with EDTA, heparin or ACD (Acid-Citrate-Dextrose) may be used to stabilize and store Pyruvate Kinase for about 15 days at 2-8°C, before to prepare the haemolyzate to test (see References 5).

• Erythrocytes 1<sup>st</sup> step: you must determine the number of erythrocytes in the normal way, reporting the result

**En = n x 10<sup>9</sup> erythrocytes/mL blood for each sample.**

• Erythrocytes 2<sup>nd</sup> step: prepare the haemolyzate. Wash 0,2 mL of whole blood with 2 mL of saline. Centrifuge 3000 rpm x 10 min. Eliminate the supernatant and resuspend with 2 mL of saline. Repeat 3x. After the elimination of the last supernatant, the washed and precipitated by centrifugation Erythrocytes must be suspended with 2,0 mL of DISTILLED WATER. Leave to wait 15 min x 2-8°C. Then centrifuge 3000 rpm x 10 min.

The supernatant is the HAEMOLYZATE. To test immediately.

A 50% of activity may be lost after 24 h at 2-8°C.

#### WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

#### ANALYTICAL PROCEDURE FOR SERUM AND URINE

##### ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against Blank Reagent
- Temperature: 37°C
- Method: end-point
- Reaction: 3 + 10 minutes
- Sample/reagent: 1/300/0.5

**Let reagents reach the working temperature before use.**

Pipette into a test tube or cuvette labelled:

BL (Blank); C (Calibrator); S: Sample (S=serum, E=erythrocytes);

	R/B	ST	S	E
<b>WORKING REAGENT</b>	2500 µL	2500 µL	2500 µL	2500 µL
Distilled water	500 µL	----	----	<b>400 µL</b>
<b>R4 - CAL</b>	----	500 µL	----	----
Sample serum	----	----	500 µL	----
Erythrocytes Sample	----	----	----	<b>100 µL</b>

Mix kindly and put in stand-by at the right temperature. Wait exactly 5 minutes; make the FIRST reading (E1 and Eb1).

Put again in stand-by at the right temperature. Wait exactly 10 minutes; make the SECOND reading (E2 and Eb2). Add:

<b>WORKING STARTER</b>	100 µL	100 µL	100 µL	100 µL
------------------------	--------	--------	--------	--------

Mix kindly and make the THIRD reading (E3 and Eb3).

Put in stand-by at the right temperature.

Wait again 10 min.; make the FOURTH reading (E4 and Eb4).

Determine the diff. of absorbance for 10 minutes due to PK activity:

**SERUM**  $\Delta A(S)$  (10 min.) = (E3 - E4) - (E1 - E2)

**ERYTHROCYTES**  $\Delta A(E)$  (10 min.) = (E3 - E4) - (E1 - E2)

**BLANK**  $\Delta A(B)$  (10 min.) = (Eb3 - Eb4) - (Eb1 - Eb2)

**CALIBRATOR**  $\Delta A(C)$  (10 min.) = (E3 - E4) - (E1 - E2)

Subtract the Blank ( $\Delta A(B)$ ) from the Calibrator and samples:

**SERUM**  $\Delta A$  (10 min.) =  $\Delta A(S)$  (10 min.) -  $\Delta A(B)$  (10 min.)

**CALIBRATOR**  $\Delta A$  (10 min.) =  $\Delta A(C)$  (10 min.) -  $\Delta A(B)$  (10 min.)

**ERYTHROCYTES**  $\Delta A$  (10 min.) =  $\Delta A(E)$  (10 min.) -  $\Delta A(B)$  (10 min.)

## ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analyzers could be completely different by what has been developed as manual determination.

## REFERENCE VALUES

Normal Values Pyruvate Kinase:

- serum 0-45 U/L (see References 6)
- erythrocytes 111-406 mU/ 10<sup>9</sup> erythrocytes /mL (see References 7)
- hemoglobin 41-152 U/g hemoglobin (see References 7)

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

## CALCULATION

Use this general formula to calculate the concentration:

### Serum

PK (U/L) =

$$\frac{(\text{Serum}) \Delta A (10 \text{ minutes})}{(\text{Calibrator}) \Delta A (10 \text{ minutes})} \times R4 - \text{CAL concentration}$$

### Erythrocytes

PK (mU/10<sup>9</sup> erythrocytes/mL = U/10<sup>9</sup> erythrocytes/L =

(Erythrocytes) ΔA(10 minutes)

$$\frac{\text{-----}}{(\text{Calibrator}) \Delta A (10 \text{ minutes})} \times R4 - \text{CAL concentration} \times 10 \times 1/n \times F$$

where 10 = dilution factor of haemolyzate  
 n = erythrocytes number mentioned in 10<sup>9</sup>  
 F = 5 [referred to the dilution 1/5 of the sample  
 (100 μl sample + 400 μl distilled water)]

But the PK in erythrocytes may be also expressed as U/g Hb. At first must be determined the concentration of Hemoglobin in g/dL. After that use this formula:

### PK (U/g Hemoglobin)

(Erythrocytes) ΔA(10 minutes)

$$\frac{\text{-----}}{(\text{Calibrator}) \Delta A (10 \text{ minutes})} \times R4 - \text{CAL conc.} \times 10 \times 1/10 \text{ Hb} \times F$$

where 10 = dilution factor of haemolyzate  
 Hb (g/dL) = Hemoglobin conc. for each sample  
 10 Hb (g/L) = 10 times Hemoglobin conc  
 F = 5 [referred to the dilution 1/5 of the sample  
 (100 μl sample + 400 μl distilled water)]

So formula becomes:

### PK (U/g Hemoglobin)

(Erythrocytes) ΔA(10 minutes)

$$\frac{\text{-----}}{(\text{Calibrator}) \Delta A (10 \text{ minutes})} \times R4 - \text{CAL conc.} \times 1/\text{Hb} \times F$$

## ANALYTICAL PERFORMANCES

### (validate on MINDRAY BS300)

The performances of the Reagent PYRUVATE KINASE have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

**Method Limitations:** are not know limitations.

**Method Linearity:** the test is linear up to 90 U/L for serum and 900

mU/10<sup>9</sup> erythrocytes /mL. For concentrations ≥ 490 U/L for serum and 900 mU/10<sup>9</sup> erythrocytes /mL, it is recommended to dilute the sample 1:10 with saline solution, test again and multiply the result x 10.

**Method Sensitivity (LoD):** the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 1.3 U/L

**Interferences:** see References 2.

Interference test criterion: recovery ± 10% of initial value. No interference found on samples with:

- total bilirubin up to 20 mg/dL;
- haemoglobin up to 450 mg/dL,
- ascorbic acid up to 37.5 mg/dL.

**Within-run Precision:** determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (U/L) ± 2s	CV %
Human 1	31.6 ± 0.7	1.2
Human 2	41.7 ± 0.8	1.1

**Run-to-run Precision:** determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (U/L) ± 2s	CV %
Human 1	31.4 ± 0.8	1.3
Human 2	41.5 ± 1.2	1.4

**Accuracy:** a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation  $y = 0.9975x - 0.103$

Correlation coefficient  $r = 0.9999$   $n = 20$

## REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.
4. Beisenherz G. et al., Z. Naturforsch. 8b, 555 (1953)
5. Beutler E., Red Cell Metabolism, Grune & Stratton, New York (1971).
6. Otto, P. et al., Klin. Wschr. 42, 75 (1964).
7. Wiesmann, U. et al., Klin. Wschr. 43, 1311 (1965)

**EDMA (EDMS) CODE 11 01 01 90 00**