ALDOLASE

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BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz,

- W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) C49-A/H56-A: Collection, Handling, Transport and Storage for Body Fluids. Quick Guide.
4. Beisenherz, G., et al., Z. Naturforsch., 8(b), 555 (1953).

Codice Ramo CND W01010104

ALDOLASE

ALDOLASE QUANTITATIVE ENZYMATIC KINETIC UV ASSAY on SERUM and PLASMA
with MINDRAY BS-300 Analyzer

500 tests (1 x 100 mL)

REF ZZ ALD100

INTENDED USE

Aldolase is an enzyme that catalyses a reverse aldol reaction. The substrate, fructose 1,6-biphosphate (F-1,6-BP) is broken down into glyceraldehyde-3-phosphate and dihydroxyacetone phosphate (DHAP). This reaction is a part of glycolysis. Aldolase helps muscle turn sugar into energy. Testing for Aldolase is done to diagnose and monitor skeletal muscle diseases.

Muscle weakness may be caused by neurologic as well as muscular problems. The measurement of Aldolase levels can help pinpoint the cause. Aldolase levels will be normal where muscle weakness is caused by neurological disease, such as poliomyelitis or multiple sclerosis; but aldolase levels will be elevated in cases of muscular disease, such as muscular dystrophy (in particular Duchenne's one).

Aldolase is also found in the liver and cardiac muscle. Damage or disease to these organs, such as chronic hepatitis or a heart attack, will also increase aldolase levels in the blood, but to a lesser degree.

Aldolase levels greater than the normal values will be also due to dermatomyositis, infectious mononucleosis, liver pancreatic or prostate cancer, myocardial infarction and polymyositis.

Also, have to be considered that strenuous exercise can temporarily increase the aldolase level in a person.

PRINCIPLE

This Reagent determines the Aldolase activity in vitro; in the presence of the substrate D-Fructose-1,6-bisphosphate and of the enzymes who manage the ancillary reactions, TIM (triosephosphate isomerase)-GDH(glycerol-3P-dehydrogenase)-LDH(lactate dehydrogenase), at the end the aldolase change NADH to NAD.

The decrease of absorbance of NADH, for oxidation to NAD, is proportional to the activity of the Aldolase in the sample.

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of Aldolase higher than 27 U/L, dilute the sample 1:10 with saline solution, repeat the determination and multiply the result by 10.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Aldolase (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

G) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

REAGENTS

Components of the kit:

R1 - BUFFER

Collidine buffer 56 mmol/L pH 7.4
F-1.6-PP 3 mmol/L

R2 - NADH

NADH 0.22 mmol/L
NaN3 < 0.1%

R3 - GDH, TIM, LDH

GDH ≥ 300 U/L
TIM ≥ 4000 U/L
LDH ≥ 500 U/L

STABILITY: the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- ALDOLASE CALIBRATOR Iyo
- NORMAL CONTROL SERUM Iyo
- PATHOLOGICAL CONTROL SERUM Iyo

REF ALDCAL MB
REF OG15500
REF OG15600

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

PREPARATION OF THE WORKING REAGENT

R1 - BUFFER: dilute the content of the bottle to 50 mL with distilled water; mix gently until dissolution.

R2 - NADH: dilute the content of the vial to 2 mL with distilled water; mix gently until dissolution.

ALDOLASE

R3 - GDH, TIM, LDH: suspension ready to use.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

STABILITY OF THE WORKING REAGENT

R1 - BUFFER diluted is stable 4 weeks at 2-8°C, in the dark.

R2 - NADH diluted is stable 4 weeks at 2-8°C, in the dark.

MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Saline solution, distilled water, Controls.

Spectrophotometer or automatic analyzer for Clinical Chemistry.

SAMPLES

- Not haemolysed fresh serum.

- Not haemolysed fresh plasma (EDTA o heparina).

Samples collection in compliance with CLSI (NCCLS)

(see References 3).

The sample can be stored at 2-8°C, up to 6 days.

(see References 1).

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

ANALYTICAL PROCEDURE ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334-365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: kinetic
- Reaction: 20 minutes
- Sample/reagent: 1/12.5

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, ST: Calibrator, S: Sample

	R/B	ST	S
R1 - BUFFER diluted	2500 µL	2500 µL	2500 µL
Saline solution	200 µL	----	----
Calibrator	----	200 µL	----
Sample	----	----	200 µL
R2 - NADH diluted	50 µL	50 µL	50 µL
R3 - GDH, TIM, LDH	10 µL	10 µL	10 µL

Mix carefully and incubate for 5 minutes at 37°C.

Read the initial absorbance of the standard (Ast1) and

sample (As1) against the reagent blank (R/B1).

After 20 minutes of incubation at 37°C, read the absorbances of

calibrator (Ast2) and sample (As2) against the reagent blank

(R/B2). Calculate:

$\Delta R/B = R/B1 - R/B2$ for the reagent blank

$\Delta AST = Ast1 - Ast2$ for the calibrator

$\Delta AS = As1 - As2$ for the sample

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analyzers could be completely different by what has been developed as manual determination.

REFERENCE VALUES (see References 1)

Normal Values Aldolase: 1.2 - 8.8 U/L

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

CALCULATION

Use this general formula to calculate the concentration:

$$\text{ALDOLASE (U/L, 37°C)} = \frac{(\Delta As - \Delta R/B)}{(\Delta Ast - \Delta R/B)} \times \text{calibrator conc.}$$

ANALYTICAL PERFORMANCES

(validate on MINDRAY BS300)

The performances of the Reagent **ALDOLASE** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 27 U/L.

For Aldolase concentrations higher than 27 U/L, it is recommended to dilute the sample 1:10 with saline solution, test again and multiply the result x 10.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 0.9 U/L.

Interferences: see References 2.

Interference test criterion: recovery \pm 10% of initial value. No interference found on samples with:

- total bilirubin up to 40 mg/dL;
- haemoglobin up to 450 mg/dL,
- lipemia [Intralipid®] up to 2000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	13.9 \pm 0.6	2.0
Human 2	18.5 \pm 1.9	5.0

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	13.1 \pm 0.8	3.2
Human 2	19.3 \pm 1.2	3.1

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

ALDOLASE

Linear regression equation $y = 1.0137x + 0.153$
Correlation coefficient $r = 0.9931$ $n = 20$

REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) C49-A/H56-A: Collection, Handling, Transport and Storage for Body Fluids. Quick Guide.
4. Beisenherz, G., et al., Z. Naturforsch., 8(b), 555 (1953).

EDMA (EDMS) CODE 11 01 01 90 00