## GLUCOSE-FRUCTOSE

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## XII. BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACCC Press (2000).
3. Tietz Textbook of Clin. Chem. and Molecular Diagnostics, W.B. Saunders Co., Philadelphia (2012), (119; 154; 581).

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Codice Ramo CND W01010211

## GLUCOSE-FRUCTOSE

For *in Vitro* Diagnostic use only**GLUCOSE and FRUCTOSE QUANTITATIVE UV ASSAY on SERUM, CEREBROSPINAL FLUID (CSF), and URINE with Analizzatore MINDRAY BS-300**

580 tests (5 x 20mL)

REF NCGF8815

### I. INTENDED USE

Glucose is the main source of energy of the living tissues; it derives from diet, from hydrolysis of hepatic glycogen and from conversion of aminoacids to glucose. Concentration and metabolism of glucose are strictly regulated by insulin, a peptide hormone produced by pancreatic B-cells that also facilitates the entry of glucose into the tissues. A decrease of insulin level or activity determines an increase of glucose in the blood. Elevated plasma concentration of glucose are found in patients with diabetes mellitus (insulin dependent or no insulin dependent). Hypoglycemia is caused by prolonged fast, administration of some drugs or congenital defects of the metabolism.

Fructose and Glucose are monosaccharides; when happens a deficiency or an absence of enzyme(s) that participates to carbohydrate metabolism, may have an accumulation of these sugars, who overflow into the urine.

Fructose is a simple monosaccharides, along with glucose and galactose, that are absorbed directly into the bloodstream during digestion. After absorption it enters the hepatic portal vein and is directed toward the liver.

The initial catabolism is in analogy with that one of glucose; it passes through fructose-1-phosphate till dihydroxyacetone phosphate and glyceraldehyde-3-phosphate who can enter the gluconeogenic pathway for glucose or glycogen synthesis as well as fatty acid and triglyceride synthesis; or to be further catabolized through the glycolitic pathway to pyruvate.

A recent study found that the elevated consumption of fructose or fructose rich foods and drinks are associated to the gout, metabolic syndrome and insulin resistance.

Deep attention is focused on total energy balance in fertility; the energy for sperm motility is derived from the metabolism of fructose carried in the seminal fluid.

### II. PRINCIPLE

In the presence of ATP, NADP and G6PDH (glucose-6P-dehydrogenase) and HK (hexokinase) the glucose produces NADPH. The intensity of the UV-colour at this wavelength is proportional to the conc. of GLUCOSE in the sample. The HK is able to phosphorylate also Fructose.

Adding PGI (phosphogluc. isomerase) new NADPH is produced from the changes of fructose-6P in glucose-6P. The new increase of the UV-colour at this wavelength is proport. to the conc. of FRUCTOSE in the sample.

### III. PRECAUTION FOR USE

1. This product has been formulated for *in vitro* diagnostic use.

2. A proportional variation of the reaction volumes does not change the result.

3. DO NOT mix Reagents from different Production lots.

4. For concentration of Glucose-Fructose  $\geq 1000$  mg/L (500 mg/L for Glucose and 500 mg/L for Fructose), dilute the sample 1:4 with distilled water, repeat the determination and multiply the result by 4.

5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

#### ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Glucose and Fructose (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) The R1 - BUFFER is supplied in surplus

G) For proteic samples may be required deproteinization (spermatic liquid for instance).

H) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

I) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

### IV. REAGENTS AND MATERIALS PROVIDED

#### Kit composition:

REF NCGF8815

**R1 - BUFFER** 3 x 70 mL

Good Buffer &gt; 20 mmol/L

**R2 - NADP, ATP** 5 x 20 mL

NADP &gt; 0.2 mmol/L

ATP &gt; 2 mmol/L

**R3 - HK, G6PDH** 1 x 2.5 mL

(READY TO USE)

HK &gt; 10 U/L

G6PDH &gt; 5 U/L

**R4 - PGI** 1 x 2.5 mL

(READY TO USE)

PGI &gt; 50 U/L

NaN3 &lt; 0.1%

### V. STORAGE AND STABILITY

The Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.  
In case of damaging of the primary container organize the waste disposal.

### VI. AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- SUBSTRATE ELEVATED CONTROL Iyo **REF OG3005**
- SUBSTRATE LOW CONTROL Iyo **REF OG3006**

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

### VII. PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - NADP, ATP** with 20 mL of **R1 - BUFFER**, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

### VIII. STABILITY OF THE WORKING REAGENT

The WORKING REAGENT is stable:

- 7 days at 2-8°C;
- 40 days at -20°C. After reconstitution FREEZE only ONE TIME, DO NOT REPEAT FREEZING.

In case of microbial contamination the preparations must be eliminated.

### IX. MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Distilled water, Controls.

Perchloric acid.

Spectrophotometer or automatic analyzer for Clinical Chemistry.

### X. SAMPLES

- Not haemolysed and not lipemic fresh serum.
- Urine.
- Cerebrospinal fluid (CSF).

For samples collection see References 3).

The sample can be stored at 2-8°C, up to 6 days. (see References 1).

### XI. SAMPLES PRETREATMENT

- Use limpid samples; If necessary, filter the samples.
- Use samples at neutral pH; neutralize the Acidic Samples at a neutral pH (max pH ≤ 8.0).

### XII. WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

### XIII. ANALYTICAL PROCEDURE ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 8 + 18 minutes
- Sample/reagent: 1/40/1/1

**Let reagents reach the working temperature before use.**

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, S: Sample

	R/B	S
<b>WORKING REAGENT</b>	1000 µL	1000 µL
Distilled water	25 µL	----
Sample	----	25 µL

Mix and incubate for about 3 minutes at 37°C. Measure the absorbance AS0 and AR/B0. For the GLUCOSE then add:

<b>R3 - HK, G6PDH</b>	25 µL	25 µL
-----------------------	-------	-------

Mix and wait the end of the react. (5-10 min.). Read AS1 and AR/B1. For the FRUCTOSE start the reactions with:

<b>R4 - PGI</b>	25 µL	25 µL
-----------------	-------	-------

Mix carefully and wait the end of the reaction (5-15 min.).

Read AS2 and AR/B2.

Calculate for the Sample Glucose

$$ASG = (AS1 - AS0);$$

Calculate for the Reagent/Blank Glucose

$$AR/BG = (AR/B1 - AR/B0).$$

**Calc. the difference for GLUCOSE  $\Delta AG = ASG - AR/BG$ .**

Calculate for the Sample Fructose

$$ASF = (AS2 - AS1);$$

calculate for the Reagent/Blank Fructose

$$AR/BF = (AR/B2 - AR/B1).$$

**Calc. the difference for FRUCTOSE  $\Delta AF = ASF - AR/BF$ .**

### ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI and MINDRAY systems. The applications on automatic analyzers could be completely different by what has been developed as manual determination.

### XIV. REFERENCE VALUES (see References 1)

Normal Values GLUCOSE adults:

- serum	740 - 1060 mg/L
- urine	< 500 mg/24h
- cerebrospinal fluid (CSF)	400 - 700 mg/L

Normal Values FRUCTOSE adults:

- serum	10 - 60 mg/L
- urine	< 30 - 65 mg/24h

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

### XV. CALCULATION FOR GLUCOSE

Use this general formula to calculate the concentration:

$$\text{Glucose conc. (g/L)} = V/v \times 1/\epsilon d \times MW/1000 \times \Delta AG$$

V = total test volume = 1.050 ml

v = sample volume = 0.025 ml

d = pathlength = 1 cm

$\epsilon$  = molar coeff. NADH = 6.3 L / mmol x cm

MW = glucose MW = 180.16

so it becomes:

$$\text{Glucose conc. (g/L)} = 1.201 \times \Delta AG$$

### XVI. CALCULATION FOR FRUCTOSE

Use this general formula to calculate the concentration:

$$\text{Fructose conc. (g/L)} = V/v \times 1/\epsilon d \times MW/1000 \times \Delta AF$$

V = total test volume = 1.075 ml

v = sample volume = 0.025 ml

d = pathlength = 1 cm

$\epsilon$  = molar coeff. NADH = 6.3 L / mmol x cm

MW = fructose MW = 180.16

so it becomes:

$$\text{Fructose conc. (g/L)} = 1.230 \times \Delta AF$$

### XVII. ANALYTICAL PERFORMANCES (Validate on MINDRAY BS300)

The performances of the Reagent **GLUCOSE-FRUCTOSE** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

**Method Limitations:** are not know limitations.

**Method Linearity:** the test is linear up to 1000 mg/L (500 mg/L for Glucose and 500 mg/L for Fructose). For concentration of Glucose-Fructose  $\geq 1000$  mg/L (500 mg/L for Glucose and 500 mg/L for Fructose), dilute the sample 1:4 with distilled water, repeat the determination and multiply the result by 4.

**Method Sensitivity (LoD) GLUCOSE:** the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 5 mg/L.

**Method Sensitivity (LoD) GLUCOSE+FRUCTOSE:** the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 13 mg/L.

**Interferences on GLUCOSE method:** see References 2.

Interference test criterion: recovery  $\pm 10\%$  of initial value. No interference found on samples with:

- total bilirubin up to 10 mg/dL;

- ascorbic acid up to 12.5 mg/dL.

**Interferences on GLUCOSE+FRUCTOSE method:** see References 2.

Interference test criterion: recovery  $\pm 10\%$  of initial value. No interference found on samples with:

- total bilirubin up to 10 mg/dL;

- ascorbic acid up to 12.5 mg/dL.

**Within-run Precision GLUCOSE:** determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (mg/L) ± 2s	CV %
Human 1	501 ± 8	0.7
Human 2	245 ± 6	1.3

**Run-to-run Precision GLUCOSE:** determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (mg/L) ± 2s	CV %
Human 1	493 ± 7	0.7
Human 2	243 ± 4	0.9

**Within-run Precision GLUCOSE+FRUCTOSE:** determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (mg/L) ± 2s	CV %
Human 1	893 ± 8	0.5
Human 2	444 ± 6	0.8

**Run-to-run Precision GLUCOSE+FRUCTOSE:** determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (mg/L) ± 2s	CV %
Human 1	887 ± 10	0.6
Human 2	445 ± 6	0.8

**Accuracy for GLUCOSE:** a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation  $y = 1.0028x - 0.007$

Correlation coefficient  $r = 0.9998$   $n = 20$


**Accuracy for FRUCTOSE:** a group of 20 sera has been tested using this procedure and using a similar reagent available on the market.

The comparison gave these results:

Linear regression equation  $y = 1.0061x - 0.014$

Correlation coefficient  $r = 0.9998$   $n = 20$

## XII. REFERENCES (see Italian version)

 IVD	In Vitro Diagnostic Medical Device	 Temperature limitation	 LOT	Batch code (XXX)	 Manufacturer	 Keep dry	 Non-sterile
 Consult Instructions for use	Use by (year/month)	 REF	Catalogue number	 Do not reuse	 Fragile, handle with care	 Keep away from heat	

EDMA (EDMS) CODE 11020111