DESTINAZIONE D'USO

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BIBLIOGRAFIA

- Anderson, N.L., et al, Cumitech 3B.; Quality Systems in the Clinical Microbiology Laboratory, Coordinating ed., A.S. Weissfeld. *American Society for Microbiology, Washington, D.C.*
- Arnold WM, Weaver RH. Quick microtechniques for the identification of cultures. *Journal of Laboratory and Clinical Medicine* 1948; 33:1334-7.
- Balzevic, D.J. and G.M. Edrer. (1975). Principles of Biochemical Tests in Diagnostic Microbiology. *John Wiley & sons, New York, NY.*
- Chapin, K. C., and T.-L. Lauderdale. 2003. Reagents, stains, and media: bacteriology, p. 354-383.
- Cruickshank, R, J. P. Duguid, B. P. Marmion, R.H.A. Swain. The Practice of Medical Microbiology. 12th Edition v2.
- Gurr, E. A Practical Manual of Medical and Biological Staining Technique. 1953
- Isenberg HD, Ed. Clinical microbiology procedures handbook, v1 Washington, DC: ASM, 1992.
- Jorgensen et al. Manual of Clinical Microbiology, *American Society for Microbiology, Washington, D.C.*
- Lennette. Manual of Clinical Microbiology. *American Society for Microbiology, Washington, D.C., 1974.*
- Lowrance, B.L., Reich, P., and Traub, W.H. (1969). *Journal of Applied Microbiology* 17:923-924.

- Murray, P. R. E. J. Baron, J. H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.), *Manual of clinical microbiology*, 8th edition. ASM Press, Washington, DC. 2003.
- Neelson, F. 1883. Ein Casuistischer Beitrag zur Lehre von der Tuberkulose. *Centralbl. Med. Wiss.* 21:497-501.
- Sutter, V.L. and W.T. Carter. (1972). *American Journal of Clinical Pathology.* 58:335-338
- Wacko, R. and J.C. Sherris. (1963). *American Journal of Clinical Pathology.* 39:429-432.

	= Use by
	= Lot number
	= Catalogue number
	= Manufacturer
	= Authorized Representative in the European Community
	= Contains sufficient for <n> tests
	= In vitro diagnostic medical device
	= Temperature limitation
	= Consult instructions for use

EC REP Advena Ltd. Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013, Malta.

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