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CALCULUS ANALYSIS III

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INSTRUCTION FOR USE

CALCULUS ANALYSIS III

READY TO USE REAGENTS FOR SEMIQUANTITATIVE COLORIMETRIC DETERMINATION OF CARBONATE, CALCIUM, OXALATE, AMMONIUM, PHOSPHATE, MAGNESIUM, URIC ACID AND CYSTINE ON KIDNEY STONES BY MANUAL METHOD

1 – INTRODUCTION AND INTENDED USE

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Kidney stones are composed of insoluble organic or inorganic substances: calcium oxalates and/or phosphates, cystine, uric ac., etc.

The origins of these conditions can be: metabolic dysfunction, obstructions, infections, hypersaturation, loss/lack of inhibitors. Almost all of the renal calculations are made up of several components, due to co-precipitation or subsequent precipitations.

The kit include all the reagents for semiquantitative colorimetric determination of Carbonate, Calcium, Oxalate, Ammonium, Phosphate, Magnesium, Uric acid and Cystine on kidney stones by manual method.

2 - PRINCIPLE OF THE METHOD

A small amount of calculus is crushed and reduced to powder, and then brought into solution. On different similar rates of this solution, they are determined semi-quantitatively:

calcium, by titration;

oxalates, ammonium, phosphates, magnesium, uric acid and cystine by colorimetric comparison.

The calculus composition is derived from the values obtained with the ruler attached to the kit. The kit contains a control sample for testing purposes.

3 - MATERIALS PROVIDED – PACKAGING

Product	Type	REF	Pack
CALCULUS ANALYSIS III	Chemical- colorimetric	3914003 (50 tests)	Various reagents contained in appropriate containers. See below for the composition of the reagents. Secondary packaging: cardboard box.

R1 – H₂SO₄

1 x 11 mL

Sulphuric Acid <100%

DANGER.

H314; P260-P264-P280-P301+P330+P331-P303+P361+P353 – P304 + P340- P305+P351+P338-P310-P321-P363-P405-P501

Note: a possible blackening to the dropping funnel doesn't influence the final result.

R2- NaOH

1 x 13 mL

Sodium hydroxide ≤ 30%

DANGER.

H314; P260-P264-P280-P301+P330+P331-P303+P361+P353 – P304 + P340- P305+P351+P338-P310-P321-P363-P405-P501

R3- Ac. Calconcarboxylic in alcohol H225

1 x 17 mL

Ac. Calconcarboxylic in alcohol H225

WARNING.

H319; P280-P337+P313

R4- EDTA

2 x 15 mL

EDTA

R5 – Borate Buffer

1 x 13 mL

Sodium hydroxide ≤2%,

Sodium Tetraborate decahydrate <5%

DANGER.

H315-H319-H360FD; P201-P202-P264-P280-P302+P352-P305+P351+P338-P308+P313-P321-P332+P313-P337+P313-P362+P363-P405-P501

R6 – FeCl₃

1 x 8 mL

Iron chloride

WARNING.

H290; P390

R7 – C₇H₆O₆S

1 x 7 mL

Sulphosalicylic Acid

WARNING.

H315-H319; P280-P337+P313

R8 – K2HgI4

Potassium Tetraiodomercurate 100%

1 x 8 mL

DANGER.

H301-H314-H341-H373-H410; P201-P202-P260-P264-P270-P273-P280-P301+P310-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

R9 – (NH4)6Mo7O24

Ammonium molybdate

1 x 13 mL

WARNING.

H315-H319; P280-P337+P313

R10 – Reducing Sol.

Methylaminophenolsulphate <1%, sodium disulphite

1 x 13 mL

WARNING

H319 – H412; P264-P273-P280-P305+P351+P338-P337+P313-P501

R11 – Buffer Sol.

Buffer
Boric acid <3%, sodium hydroxide <1%

1 x 25 mL

DANGER.

H315-H319-H360FD;P201-P202-P264-P280-P302+P352-P305+P351+P338-P308+P313-P362+P364-P405-P501

R12 – Dye Sol.

Blue dyeing agent

1 x 25 mL

WARNING.

H319; P280-P337+P313

R13 – Phosphomolibdic

Phosphomolibdic <20%

1 x 8 mL

DANGER.

H314;P260-P264-P280-P301+P330+P331-P303+P351+P338-P310-P321-P363-P405-P501

R14- NH4OH

Ammonium hydroxide <50%

1 x 25 mL

DANGER.

H314–H335-H410;P260-P261-P264-P271-P273-P280-P301+P330+P331–P303+P361+P353–P304+P340-P305+P351+P338–P310–P312-P321-P363-P391-P403+P233-P405-P501

R15 – Reducing Powder

Na2SO3

1 x 25 g

R16- Nitroprussiate

Sodium nitroprussiate

1 x 25 g

RSC- CONTROL

CONTROL for semiquantitative determination. Expected values on attached sheet.

1 x 1 g

WARNING.

H319; P280-P337+P313

CR – Calculator Rule

N°1

CS - Chromatic Scale

N°1

TUBES 15 - Sample Tubes 15 mL

N°5

TUBES 5 – Test Tubes 5mL

N°7

TUBE 50 – Test Tube 50 mL

N°1

DOS 25 – Dispenser O.25

N°2

DOS 20 – Dispenser O.20

N°2

Note: All the plastics used, from TUBES 15 to DOS 20 can be washed with distilled water, placed to dry completely to be reused.

HAZARD STATEMENTS

H290 May be corrosive to metals.

H301 Toxic if swallowed.

- H314** Causes severe skin burns and eye damage.
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.
H360FD May damage fertility. May damage the unborn child.
H341 Suspected of causing genetic defects.
H373 May cause damage to organs through prolonged or repeated exposure.
H410 Very toxic to aquatic life with long lasting effects.
H412 Harmful to aquatic life with long lasting effects.

PRECAUTIONERY STATEMENTS

- P201** Obtain special instructions before use.
P202 Do not handle until all safety precautions have been read and understood.
P260 Do not breathe dust/fume/gas/mist/vapours/ spray.
P261 Avoid breathing dust/fume/gas/mist/vapours/ spray.
P264 Wash ... thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
P271 Use only outdoors or in a well-ventilated area.
P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection/ face protection.
P301+P310 IF SWALLOWED: Immediately call a POISON CENTER/doctor/...
P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P302+P352 IF ON SKIN: Wash with plenty of water/...
P303+P361+P353 If on skin (or hair): remove/take off immediately all contaminated clothing. Rinse skin with water/shower.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313
P310 Immediately call a poison center or doctor/physician.
P312 Call a POISON CENTER/doctor/.../if you feel unwell.
P314 Get medical advice/attention if you feel unwell.
P321 Specific treatment (see ... on this label).
P330 Rinse mouth.
P332+P313 If skin irritation occurs: Get medical advice/attention.
P337+P313 If eye irritation persists: Get medical advice/ attention.
P362+P364 Take off contaminated clothing and wash it before reuse.
P363 Wash contaminated clothing before reuse.
P390 Absorb spillage to prevent material damage.
P391 Collect spillage.
P403+P233 Store in a well-ventilated place. Keep container tightly closed.
P405 Store locked up.
P501 Dispose of contents/container to ...

4 - MATERIALS REQUIRED BUT NOT PROVIDED

Clean porcelain pestle and mortar. Spatula. Distilled water.

5 - PRECAUTIONS AND WARNINGS

- This product has been formulated for in vitro diagnostic use.
- DO NOT mix Reagents from different Production lots.
- It is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.
- Do not use after expiry date. The reagents' quality cannot be guaranteed beyond their shelf-life date or if the reagents are stored under inappropriate conditions.

ATTENTION!

- The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.
- The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.
- A lot of factors, as room temperature, the working reagent temperature, may affect the tests performances.
- All the precautions normally used in the laboratory must be respected for reagents handling.
- Follow exactly what described here to have right results.
- The Certificates of Analysis and the Safety Data Sheet of the product are available on the website. www.masciabrunelli.it.
- The information provided in this document has been defined to the best of our knowledge and ability and represents a guideline for the proper use of the product but without obligation or liability. In all cases existing local laws, regulations and standard procedures must be observed for the examination of samples collected from human and animal organic districts, for environmental samples and for products intended for human or animal consumption. Our information does not relieve our customers from their responsibility for checking the suitability of our product for the intended purpose.
- Notify Mascia Brunelli Spa and the Relevant Authorities of any serious incidents occurring in connection with the in vitro diagnostic device. complaint@masciabrunelli.it

6 - STORAGE CONDITIONS AND SHELF LIFE

The Reagents are stable up to the expiry date mentioned on the labels, stored at 15-25°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

STABILITY AFTER FIRST OPENING: The product is stable up to the expiry date mentioned on the labels after the first open if stored at 15-25°C.

7 – QUALITY CONTROL

Using the control in the kit, as it contains all analytes that can be determined with this device; you can verify the chemical reaction for your reference.

8 – PREPARATION OF THE WORKING REAGENT

Ready-to-use. Mix kindly before use and let the reagent reach the working temperature before use. Close immediately after handling. The Reagents have to be used correctly, to avoid contamination.

During the time and/or due to temperature changes, the R8 and R13 may form a suspension or a precipitate; this does not affect the final result of the analysis.

9 – SAMPLES

Renal or urinary stone.

10 – WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

11 – ANALYTICAL PROCEDURE

PRELIMINARY PHASE AND TEST FOR CARBONATE

After an accurate organoleptic examination (form, colour, etc.), crush the kidney stone and reduce to powder in a mortar. With a dosing spoon (ref. DOS20), transfer a sample of the well mixed powder (a level scoop, 15-20mg) in a Sample tube (ref. TUBES15).

Use the ref. RSC – CONTROL as a sample.

Add 5 drops of R1 - H₂SO₄ and close the cap. Shake gently by holding it upright or tapping it on the bottom with your fingers to mix; do not tip it upside down to avoid losing material. Place in a tube holder until completely dissolved. **The formation of effervescence during dissolution and/or agitation, indicates the presence of carbonates.**

Add about 40 mL of distilled water in the TUBE 50. Unscrew the cap of TUBES 15 with the dissolved sample and adding about 10-15mL of distilled water from TUBE 50; close and mix 2-3 times by inversion. Then mix the contents into the TUBE 50. Repeat the procedure once again to make sure that you have brought all the dissolved sample into the TUBE 50.

Then make up to 50mL with distilled water and screw in the cap; mix for inversion 2-3 times.

Distribute the contents of the TUBE 50 in the 7 available TUBES 5, at a level of 5mL for 6 tubes (5mL mark) and placing only 1mL exactly in the 7th tube (for Magnesium). That will then be taken to the 5mL mark with distilled water.

CALCIUM DETERMINATION

Principle

Calcium is determined in a complexometric way with EDTA, with calconcarboxylic acid as an indicator.

Procedure

Add to a tube with sample to the 5 mL mark

- R2 2 drops and

- R3 3 drops (4 if you want a more intense colour) then add

- R4 drop by drop, until the colour of the solution changes from red to blue. Count exactly the number of drops added.

The Calcium percentage in Calcium is obtained by multiplying the number of drops x 5.

OXALATE DETERMINATION

Principle

Oxalates discolorize the coloured complex formed by sulphosalicylic acid and Fe (III).

Procedure

Add to a tube with sample to the 5 mL mark, closing and shaking gently in succession

- R5 2 drops,

- R6 3 drops and

- R7 3 drops; close and leave to stand for 2 minutes.

Classify by comparing the colour of the solution in the reaction tube with its colorimetric scale (CS). Always look through the solution from top to bottom and read the percentage of oxalate in the calculation. Assess intermediate grades as well.

AMMONIUM DETERMINATION

Principle

Sodium tetraiodomercurate and ammonium ion form yellow to brown coloured solutions.

Procedure

Add to a tube with sample to the 5 mL mark, closing and shaking gently in succession

- R8 3 drops,

- R2 3 drops.

Classify by comparing the colour of the solution in the reaction tube with its colorimetric scale (CS). Always look through the solution from top to bottom and read the percentage of ammonium in the calculation. Assess intermediate grades as well.

PHOSPHATES DETERMINATION

Principle

The phosphomolybdic acid formed by the addition of ammonium molybdate is reduced with a suitable reducing agent to molybdenum blue.

Procedure

Add to a tube with sample to the 5 mL mark, closing and shaking gently in succession

- R9 5 drops,

- R10 5 drops, then leave to rest for 5 minutes.

Classify by comparing the colour of the solution in the reaction tube with its colorimetric scale (CS). Always look through the solution from top to bottom and read the percentage of phosphates in the calculation. Assess intermediate grades as well.

MAGNESIUM DETERMINATION

Principle

In buffered solution magnesium reacts with 1-azo-2-hydroxy-3-(2,4-dimethylcarboxanilide)-naphthalene-1'-(2-hydroxybenzene-5-sulfonate sodium), forming a red complex.

Procedure

In the **TUBESS** containing 1 mL of sample, make up to the 5 mL mark with distilled water, tap and mix. Add closing and shaking gently in succession

- **R11** 10 drops,
- **R12** 10 drops. Wait 1 minute.

Classify by comparing the colour of the solution in the reaction tube with its colorimetric scale (**CS**). Always look through the solution from top to bottom and read the percentage of magnesium in the calculation. Assess intermediate grades as well.

URIC ACID DETERMINATION

Principle

In buffered solution, uric acid reduces phosphomolybdenic acid to molybdenum blue.

Procedure

Add to a tube with sample to the 5 mL mark

- **R13** 3 drops, waiting 2 minutes.

Then add

- **R5** 2 drops and mixing.

Classify immediately by comparing the colour of the solution in the reaction tube with its colorimetric scale (**CS**). Always look through the solution from top to bottom and read the percentage of uric acid in the calculation. Assess intermediate grades as well.

This comparison must be made within 10" of the addition of the R5 reagent, as the staining is not stable and tends to turn blue.

CYSTINE DETERMINATION

Principle

Sodium sulfite reduces cystine to cysteine, which in alkaline solution with sodium nitroprussiate forms a red colouring.

Procedure

Add to a tube with sample to the 5 mL mark, closing and shaking gently in succession

- **R14** 10 drops,
- **R15** 1 spoonful with **DOS 25** well levelled of powder, stirring until completely dissolved.

After 1' add

- **R16** 1 spoonful with **DOS 25** well levelled of powder, stirring until completely dissolved.

After 30" classify by comparing the colour of the solution in the reaction tube with its colorimetric scale (**CS**). Always look through the solution from top to bottom and read the percentage of cystine in the calculation. Assess intermediate grades as well.

NOTE

It is possible to reduce the time of analysis by performing the determination of ammonium in the 2' waiting time necessary for the determination of oxalate and magnesium during the 5' waiting time necessary for the determination of phosphates.

12 – CALCULATION

The chemical elements determined with this kit generally form the following compounds in kidney stones, postulated as the basis of the calculator rule (**CR**):

$\text{CaC}_2\text{O}_4 \times \text{H}_2\text{O}$	Whewellite	Calcium oxalate
$\text{MgNH}_4\text{PO}_4 \times 6 \text{H}_2\text{O}$	Struvite	Magnesium ammonium phosphate
$\text{CaHPO}_4 \times 2 \text{H}_2\text{O}$	Brushite	Dicalcium phosphate
$\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$	Apatite	Tricalcium phosphate
Ammonium urate		
Uric acid		
Cystine		

Weddellite ($\text{CaC}_2\text{O}_4 \times 2 \text{H}_2\text{O}$) and Apatite carbonate ($\text{Ca}_{10}[\text{PO}_4]_6\text{CO}_3$) have been overlooked because such differentiation has no clinical significance.

From the determined percentages of the single components we refer to the composition of the kidney stone as follows:

- for each cation, the possible compounds of the above mentioned list are hypothesized, based on the anions really present;
- Using the calculator rule (**CR**) you can then obtain the quantities of the single compounds probably present.

1) Calcium Oxalate (Whewellite):

- 1.1 set on the scale of oxalates the percentage of oxalate determined experimentally and read the corresponding value on the scale of the calcium oxalate;
- 1.2 trace the corresponding amount of calcium;
- 1.3 if the amount of calcium determined experimentally is greater than the quantity calculated as described above, subtract the amount of non-combined calcium as oxalate (see Calcium Phosphate).

2) Magnesium ammonium phosphate (Struvite):

- 2.1 set on the scale of magnesium the percentage of magnesium determined experimentally and read the corresponding value on the scale of the Struvite corresponding to magnesium ammonium phosphate;
- 2.2 trace the corresponding quantities of ammonium and phosphate;
- 2.3 if the amount of ammonium or phosphate determined experimentally is greater than the quantity calculated as described above, subtract the amount of ammonium and noncombined phosphate as magnesium ammonium phosphate (see ammonium urate for ammonium, calcium phosphates for phosphate).

3) Ammonium Urate

- 3.1 set on the scale of ammonium the percentage of ammonium determined experimentally, or obtained by subtraction (see point 2.2) and read the corresponding value on the ammonium urate scale;
- 3.2 trace the corresponding amount of uric acid;
- 3.3 if the amount of experimentally determined uric acid is greater than the amount calculated as described above, subtract the residual amount of uric acid.

4) Calcium phosphate

- 4.1 set on the scale of calcium the percentage of calcium determined experimentally, or that obtained by subtraction (see point 1.2);
- 4.2 check on the scale of phosphates which are the values obtained by calculation who correspond better to the values determined experimentally or subtracted, as described in point 2.2;
- 4.3 Finally, trace the brushite concentration from the respective scales.

Example

From stone analysis we obtained the following concentrations:

calcium 35%, oxalate 15%, phosphate 40%
 Kidney stone may therefore contain calcium oxalate and calcium phosphates.
 15% of oxalate correspond to 25% of calcium oxalate on calculator rule (CR).
 25% of calcium oxalate correspond to 7% of calcium.

Not combined as oxalate therefore remains 28% calcium.
 28% calcium corresponds to 40% phosphates on calculator rule (CR): 70% apatite is thus deducted from the lower scale of calculator rule (CR).

It is concluded that the stone consists of 28% calcium oxalate and 70% apatite

Note: For method reasons, the sum of components does not generally give a value exactly equal to 100%.

13 – COMPARISON OF METHODS

The comparison of the results obtained with CALCULUS ANALYSIS III and a similar reagent available on the market has given almost identical results.

14 – BIBLIOGRAPHY

1. Shokouhi E. et al., Res.J.Biol.Sci. 3, 620 (2008).
2. Iqbal M.W. et al., CJPAS 2, 139 (2008).
3. Lopez M. et al., Pediatr. Nephrol. 25, 49 (2010).
4. Pandeya A. et al., Nepal Med. Coll. J. 12, 190 (2010).
5. Sandhya A. et al., IJABPT 1, 175 (2010).
6. Raina A.F. et al., Int.J.Adv.Med. 4, 1477 (2017).
7. Cho S.T. et al., Int. J. Urology 20, 208 (2013)

TABLE OF APPLICABLE SYMBOLS

	In Vitro Diagnostic Medical Device		Temperature limitation		Batch code (XXX)		Manufacturer		Keep dry		Unique device identifier
	Consult Instructions for use		Use by (year/month)		Catalogue number		Do not reuse		Fragile, handle with care		Keep away from heat

REVISION HISTORY

Version	Description of changes	Date
Instructions for Use (IFU) - Revision 2	Updated layout, content and Symbols in the SDS	2024/03

Note: minor typographical, grammatical, and formatting changes are not included in the revision history.