## AMMONIA

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## BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) C49-A/H56-A: Collection, Handling, Transport and Storage for Body Fluids. Quick Guide.

**Codice Ramo CND W01010301**

## AMMONIA

QUANTITATIVE KINETIC UV DETERMINATION OF PYRUVATE KINASE (PK)  
in SERUM and ERYTHROCYTES on MINDRAY BS-300 Analyzer

300 tests (20 x 3 mL)

**REF** NAAM8870

### INTENDED USE

Deamination is a fundamental step in the proteic metabolism. During it, the aminoacids lose the amino group (NH<sub>2</sub>). This group pass to alpha-ketoglutarate (Krebs cycle intermediate), forming glutamate; this last reacts by oxidative deamination in the mitochondrial matrix giving free ammonia (NH<sub>3</sub>).

Ammonia is a toxic molecule, particularly for brain.

An healthy organism is able to change it in atoxic compound, which are also transport and pre-elimination means:

- 1) may be add to alpha-ketoglutarate to form again glutamate,
- 2) may be incorporated in a glutamate molecule to give glutamine,
- 3) may be used for the carbamylphosphate synthesis, urea cycle precursor, which transforms ammonia into urea at hepatic level; urea is no toxic molecule and is then eliminated by urine.

In presence of serious hepatic diseases such as cirrhosis, hepatitis and Reye's Syndrome, the Ammonia is accumulated in the blood. Elevated ammonia concentrations are toxic, especially for the brain and can be involved in the development of the portosystemic encephalopathy.

### PRINCIPLE

Ammonia reacts with 2-glutarate to give L-glutamate, in the presence of glutamate dehydrogenase (GLDH) and reduced nicotinamide-adenine dinucleotide phosphate (NADPH); in this reaction NADPH is oxidized.

The amount of NADPH oxidized at UV wavelength is proportional to the concentration of Ammonia in the tested sample.

Using the standard contained in the kit it is possible to prepare a Calibration Curve to refer.

Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample.

### PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of Ammonia higher than the 1000 µg/dL (587.0 µmol/L), dilute the sample 1:2 with saline solution, repeat the determination and multiply the result by 2.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the

material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

### ATTENTION!

A) Applications on routine analysers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyser.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Ammonia (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) To work in a free atmosphere from ammonium vapors (DO NOT SMOKE in laboratory).

G) NEVER FORGET that Ammonia evaporates very easily from samples.

H) Keep always the samples in test-tube well closed and in ice.

I) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

J) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

### REAGENTS

Components of the kit:

**R1 - BUFFER** **REF** AM8870 **1 x 60 mL**

Good buffer > 10 mmol/L

**R2 - ketoglutarate** **20 x 3 mL**

ketoglutarate > 4 mmol/L

**R3 - GLDH, NADPH** **1 x 16 mL**

(READY TO USE)

GLDH > 200 U/L ADP > 3 mmol/L

NADPH > 0.1 mmol/L Eccipienti

**R4 - CAL** **1 x 5 mL**

(READY TO USE)

Solution of Ammonia = 500 µg/dL (= 293.5 mmol/L)

NaN<sub>3</sub> < 0.1%

**STABILITY:** the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8 °C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

## QUALITY CONTROL

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

## PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - ketoglutarate** with 3 mL of **R1 - BUFFER**, mix gently avoid foaming.

**R3 - GLDH, NADPH** is ready to use. Before use, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

## STABILITY OF THE WORKING REAGENT

The WORKING REAGENT is stable 10 days at 2-8°C. Till 40 days frozen at -20°C. FREEZE only ONE TIME. DO NOT REPEAT FREEZING.

## MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 500 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Saline solution.

Spectrophotometer or automatic analyser for Clinical Chemistry.

## SAMPLES

• NOT HAEMOLYSED fresh plasma.

Samples collection in compliance with CLSI (NCCLS)

(see References 3).

On recommend to keep samples on ICE and to separate the plasma from blood without delay (see References 1).

## WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

## ANALYTICAL SPECTROPHOTOMETER PROCEDURE ON MANUAL

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: fixed-time
- Reaction: 7 - 10 minutes
- Sample/reagent: 1/10/2.5

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, ST: Calibrator, S: Sample

	R/B	S	ST
<b>WORKING REAGENT</b>	1000 µL	1000 µL	1000 µL
Distilled water	100 µL	----	----
<b>R4 - CAL</b>	----	100 µL	----
Sample	----	----	100 µL

Mix and incubate for about 3 minutes at 37°C. Then add:

<b>R3 - GLDH, NADPH</b>	250 µL	250 µL	250 µL
-------------------------	--------	--------	--------

Mix carefully, incubate at 37°C.

After 2 minutes read AS1 (sample), AST1 (standard) and AR/B1 (blank).

Wait again 5 minutes and read AS2, AST2 and AR/B2.

Calculate for the sample AS = (AS1 – AS2).

Calculate for the standard AST = (AST1 – AST2).

Calculate for the Reagent/Blank AR/B = (AR/B1 - AR/B2).

**Calculate the difference**  $\Delta AS = AS - AR/B$  and

$\Delta AST = AST - AR/B$ .

## ATTENTION!

**The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analysers could be completely different by what has been developed as manual determination.**

## REFERENCE VALUES (see References 1)

Normal Values Ammonia : 15 - 45 µg/dL (8.8 – 26.4 µmol/L).

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

## CALCULATION

Use this general formula to calculate the concentration:

$[\Delta AS / \Delta AST] \times 500^* = \text{AMMONIA } (\mu\text{g/dL})$

\*500 in R4 - **CAL** concentration

Ammonia in µmol/L = 0.587 x µg/dL of Ammonia.

## ANALYTICAL PERFORMANCES

### (validate on MINDRAY BS300)

The performances of the Reagent **Ammonia** have been tested with a MINDRAY BS300 analyser. The data, while representing the characteristics of the product, could be different for each laboratory and for different analysers.

**Method Limitations:** are not know limitations.

**Method Linearity:** the test is linear up to 1000 µg/dL

(587.0 µmol/L). For concentrations ≥ 1000 µg/dL, it is recommended to dilute the sample 1:2 with saline solution, test again and multiply the result x 2.

**Method Sensitivity (LoD):** the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 2.8 µg/dL.

**Interferences:** see References 2.

Interference test criterion: recovery ± 15% of initial value. No interference found on samples with:

- total bilirubin up to 20 mg/dL;
- lipemia [Intralipid®] up to 1000 mg/dL;
- ascorbic acid up to 50 mg/dL.

**Within-run Precision:** determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (µg/dL) ± 2s	CV %
Human 1	490.7 ± 20.1	2.1
Human 2	67.8 ± 5.1	3.7

**Run-to-run Precision:** determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (µg/dL) ± 2s	CV %
Human 1	494.3 ± 15.9	1.6
Human 2	68.2 ± 7.8	5.8

**Accuracy:** a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation  $y = 0.9988x + 2$   
Correlation coefficient  $r = 0.9997$   $n = 20$

## REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) C49-A/H56-A: Collection, Handling, Transport and Storage for Body Fluids. Quick Guide.

**EDMA (EDMS) CODE 11 03 01 01 00**