L-LACTIC ACID

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Regressione lineare $y = 0.9904x + 0.007$

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XIV. BIBLIOGRAFIA

- Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
- Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
- CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.

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L-LACTIC ACID

For *in Vitro* Diagnostic use only

L-LACTIC ACID QUANTITATIVE UV ASSAY on PLASMA, URINE and LIQUOR with MINDRAY BS-300

 Σ (488 tests (5 x 20 mL))

REF NALA8810

I. INTENDED USE

L-lactic acid, and its salt lactate, are carboxylic acid that plays a role in various biochemical processes. L-lactate is constantly produced from pyruvate via lactate dehydrogenase (LDH) during normal metabolism and exercise. It does not increase in concentration until the rate of lactate production exceeds the rate of lactate removal, which is governed by a number of factors, including monocarboxylate transporters, concentration and isoform of LDH, and oxidative capacity of tissues.

The increased lactate produced can be removed in a number of ways, including:

-oxidation to pyruvate by well-oxygenated muscle cells; it is then directly used to fuel the Krebs cycle.

-conversion to glucose via gluconeogenesis in the liver and release back into the circulation. If not released, the glucose can be used to build up the liver's glycogen stores if they are empty.

Lactate, and not glucose, is preferentially metabolized by neurons in the brain; the glial cells are responsible for transforming glucose into lactate, and for providing lactate to the neurons.

The role of lactate for brain metabolism seems to be even more important at early stages of development (prenatal and early postnatal), with lactate at these stages having higher concentrations in body liquids.

Pathological lactate production occurs when there is inadequate tissue oxygen supply, or if there is defective pyruvate clearance or increased pyruvate production. Type A lactic acidosis occurs when there is obvious tissue hypoxia and Type B when tissue hypoxia is not apparent. Hypoxia is by far the most common cause of lactic acidosis and can occur with shock from any cause. Other factors which may cause tissue hypoxia include severe asthma, severe anaemia and carbon monoxide or cyanide poisoning.

II. PRINCIPLE

The L-Lactic acid is changed in pyruvate by L-LDH (lactate dehydrogenase) in the presence of NAD. This reaction is helped by a secondary one, who transforms the pyruvate originated before in the presence of L-glutamate and GPT (glutamate-pyruvate transaminase).

The intensity of the UV-colour at this wavelength is proportional to the concentration of L-Lactic acid in the tested sample.

III. PRECAUTION FOR USE

1. This product has been formulated for *in vitro* diagnostic use.

2. A proportional variation of the reaction volumes does not change the result.

3. DO NOT mix Reagents from different Production lots.

4. For concentration of L-Lactic Acid higher than 400 mg/L (4.44 mmol/L), dilute the sample 1:4 with saline solution, repeat the determination and multiply the result by 4.

5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of L-Lactic Acid (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) The R1 - BUFFER is supplied in surplus

G) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

H) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

IV. REAGENTS AND MATERIALS PROVIDED

Components of the kit:

	NALA8810
R1 - BUFFER	3 x 70 mL
Good buffer	> 10 mmol/L
R2 - NAD, LDH	5 x 20 mL
Glutamic acid	> 3 mmol/L
NAD	> 0.4 mmol/L
LDH	> 300 U/L
R3 - GPT (ALAT)	1 x 2.5 mL
(READY TO USE)	
GPT(ALAT)	> 20 U/L
NaN ₃	< 0.1%

MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 μ L.

Disposable micropipettes tips.

Transparent glass tubes for sample dilution.

Calibrators and Controls

Spectrophotometer or automatic analyser for Clinical Chemistry.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- SUBSTRATE ELEVATED CONTROL Iyo

REF OG3005

- SUBSTRATE LOW CONTROL Iyo

REF OG3006

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

V. STORAGE AND STABILITY

the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

VI. PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - NAD, LDH** with 20 mL of **R1 - BUFFER**, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

VII. STABILITY OF THE WORKING REAGENT

The working reagent is stable 1 day at 2-8°C.

VIII. SAMPLES

- Not haemolysed and not lipemic fresh plasma, ONLY with NaF;
- Urine;
- Liquor (CSF).

Use limpid samples; If necessary, filter the samples.

Use samples at neutral pH; neutralize the Acidic Samples at a neutral pH (max pH ≤ 8.0).

Samples collection in compliance with CLSI (NCCLS) (see References 3).

The sample can be stored at 2-8°C, up to 6 days. (see References 1).

IX. WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

X. ANALYTICAL PROCEDURE ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end point
- Reaction: 8 - 16 minutes
- Sample/reagent: 1/40/1

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled: R/B: Reagent Blank, S: Sample

	R/B	S
WORKING REAGENT	1000 µL	1000 µL
Distilled water	25 µL	----
Sample	----	25 µL

Mix and incubate for about 3 minutes at 37°C. Measure the absorbance AS1, and AR/B1. Then add:

R3 - GPT (ALAT)	25 µL	25 µL
------------------------	-------	-------

Mix carefully, incubate at 37°C and wait the end of the reaction (5-10 min.). Read AS2, and AR/B2.

Calculate for the sample $AS = (AS2 - AS1)$;

Calculate for the Reagent/Blank $AR/B = (AR/B2 - AR/B1)$.

Calculate the difference $\Delta AS = AS - AR/B$.

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analysers could be completely different by what has been developed as manual determination.

XI. REFERENCE VALUES FOR WHOLE BLOOD (see References 1)

Normal Values L-Lactic acid:

Plasma Venous 4.5-19.8 mg/dL (0.5-2.2 mmol/L)

Arterial 4.5-14.4 mg/dL (0.5-1.6 mmol/L)

Urine 496-1982 mg/24 h (5,5 - 22 mmol/24 h)

CFS Adult 10 - 22 mg/dL (1,1 - 2,4 mmol/L)

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

XII. CALCULATION

Use this general formula to calculate the concentration:

$$\text{L-Lactic conc. (g/L)} = V/v \times 1/\epsilon d \times MW/1000 \times \Delta A$$

V = total test volume = 1.050 mL

v = sample volume = 0.025 mL

d = pathlength = 1 cm

ϵ = molar coeff. NADH = 6.3 L / mmol x cm

MW = lactic acid MW = 90.1

so it becomes:

$$\text{L-Lactic acid conc. (g/L)} = 0.601 \times \Delta A$$

XIII. ANALYTICAL PERFORMANCES (Validate on MINDRAY BS300)

The performances of the Reagent **L-LACTIC ACID** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 400 mg/L

(2.44 mmol/L). For concentrations \geq 400 mg/L (2.44 mmol/L), it is recommended to dilute the sample 1:4 with saline solution, test again and multiply the result x 4.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 2 mg/L (0.02 mmol/L).

Interferences: see References 2.

Interference test criterion: recovery \pm 10% of initial value. No interference found on samples with:

- total bilirubin up to 40 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	388 \pm 4	0.6
Human 2	246 \pm 2	0.5

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:


Sample	Mean (U/L) \pm 2s	CV %
Human 1	389 \pm 5	0.7
Human 2	245 \pm 5	1.1

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation $y = 0.9904x + 0.007$

Correlation coefficient $r = 0.9987$ $n = 20$

XIV. REFERENCES (see Italian version)

	In Vitro Diagnostic Medical Device		Temperature limitation		Batch code (Axxx)		Manufacturer
	Consult instructions for use		Use by (year/month)		Catalogue number		Do not reuse
	Keep dry		Non-sterile		Fragile, handle with care		Keep away from heat

EDMA (EDMS) CODE 1102011800