RAME

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RAME

For *in Vitro* Diagnostic use only

COPPER QUANTITATIVE COLORIMETRIC ASSAY on SERUM, PLASMA and URINE with MINDRAY BS-300

 Σ (2 x 50 mL) 330 tests

REF NB12000

I. INTENDED USE

Copper is present in the body in a particular area: it is a trace element associated with a large number of metalloproteins.

In the biological systems Copper is in the 2 valence states and permits to metalloproteins to play a central role in the oxidation-reduction reactions. Ceruloplasmin, ascorbate oxidase, tyrosinase, lysyl oxidase, superoxide dismutase, cytochrome c oxidase and so on, are a few of the most important copper-containing enzymes: they can bind and react directly with molecular oxygen.

Several pathologies and alteration of the physiological conditions are due to modifications or to the loss of copper-enzymes: from connective tissue defects (cardiac, skeletal and vascular) to failure in pigmentation, from ataxia in motor neurons to abnormal catecholamine conversions. After these, Copper play an important role in the Iron metabolism. Low levels may interfere in the right iron absorption and severe copper deficiency gives elevated anemia: ceruloplasmin has a strong action on ferrous iron to give ferric iron, exactly before its binding to plasma transferrin.

II. PRINCIPLE

Copper dissociated from proteins, in particular conditions of ionic strength, gives with Di-Br-PAESA a stable colored complex, which intensity of color is proportional at the concentration of Copper in the sample.

III. PRECAUTION FOR USE

1. This product has been formulated for *in vitro* diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of copper higher than 600 μ g/dL, dilute the sample 1:2 with saline solution, repeat the determination and multiply the result by 2.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However, it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

- A) Applications on routine analysers may be totally different from what developed as manual determination; in addition, the procedures are specific for each analyser.
- B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Copper (see References 2).
- C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.
- D) The clinical diagnosis cannot be done correctly using the result of only one test but have to be done integrating critically the results of different laboratory tests and clinical data.
- E) To avoid interference or contaminations, use plastic material throwaway or very clean tubes washed with diluted HCl and distilled water.
- F) The presence of a precipitate in the R1 - BUFFER over time is physiological and does not lead to degradation in the functioning and stability of the reagent.
- G) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.
- H) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.
- I) All the precautions normally used in the laboratory must be respected for reagents handling.
- All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

IV. REAGENTS AND MATERIALS PROVIDED

Kit composition:

R1 - BUFFER

Acetate Buffer 0.1 M pH 4.9; reducing agents and preservatives

R2 - DiBr-PAESA

DiBr-PAESA >0.001 mmol/L, Reducing agents, Stabilizers

R3 - CAL Standard 200 serum/plasma

(READY TO USE) Solution of Copper = 200 μ g/dL

R4 - Blank reagent urine and preservatives

(READY TO USE) Solution of Copper = 0 μ g/dL

R5 - CAL Standard 10 urine and preservatives

(READY TO USE) Solution of Copper = 10 μ g/dLREF NB12000
1 x 50 mL

1 x 50 mL

1 x 5 mL

1 x 15 mL

1 x 5 mL

MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 μ L.

Disposable micropipettes tips.

Transparent glass tubes for sample dilution.

Saline solution, Controls

Spectrophotometer or automatic analyser for Clinical Chemistry.

V. AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use):

- NORMAL CONTROL SERUM Iyo

REF NBCNU

- PATHOLOGICAL CONTROL SERUM Iyo

REF NBCPU

- COPPER CONTROLS normal/elevated

REF NBCOPC067

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

VI. STORAGE AND STABILITY

The Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations. **Component R1 - BUFFER can be stored at 15°C-25°C; if stored at 2°C-8°C, it can form precipitates. In this case, place Reagent R1 in a water bath for 5 minutes at 30-35°C and mix by inversion.** In case of damaging of the primary container organize the waste disposal.

VII. PREPARATION OF THE WORKING REAGENT

Prepare a working reagent by mixing Reagent 1 in equal parts with Reagent 2 or use the two reagents in equal parts directly in the reaction cuvette. If prepared, the working reagent is stable for 20 days at room temperature. Mix kindly before use and let the reagent reach the room temperature before use. Close immediately after handling. The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

VIII. STABILITY AFTER FIRST OPEN

The product is stable 4 months after the first open if stored at 2-8°C. Component R1 can be stored at 15 °C-25 ° C for 4 months.

IX. SAMPLES

- Not haemolysed fresh serum.
- Not haemolysed fresh plasma in heparin.
- Fresh urine, clear; in the presence of turbidity or precipitate, filter through filter paper or centrifuge for 10/15 minutes at 3000 rpm. All the urine have to be adjusted until pH 7.00.

At the calculation, remember the dilution ratio of the sample, if made. Samples collection in compliance with CLSI (NCCLS) (see References 3).

X. WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

XI. ANALYTICAL PROCEDURE ON MANUAL SPECTROPHOTOMETER

- Wavelength: 580 nm (570 - 600 nm)
- Pathlength: 1 cm
- Reading: against blank reagent
- Temperature: 37°C
- Method: end-point
- Reaction: 10 minutes
- Ratio Sample/reagent serum/plasma: 1/15
- Ratio Sample/reagent urine: 1/3.3

Let reagents reach the working temperature before use.

Pipette in vial or cuvette so labelled: R/B: Reagent Blank, ST: Standard, S: Sample.

If Reagent 1 + Reagent 2 is used, mix gently.

Then add the water, standard, sample, respectively. The table below is an example:

SERUM/PLASMA	R/B	ST _{Serum/plasma}	S
Working reagent (R1+R2)	1 (0.5+0.5) mL	1 (0.5+0.5) mL	1 (0.5+0.5) mL
Distilled water	66 µL	---	---
R3 – CAL serum/plasma	---	66 µL	---
Sample	---	---	66 µL

Mix carefully. Read the absorbance of the standard (Ast) and sample (As) against reagent blank, after 10 minutes at 37°C. The colour is stable 30 minutes at room temperature.

CALCULATION SERUM/PLASMA

$(As / Ast) \times 200 * = \mu\text{g/dL}$ of Copper

* is the **R3 - CAL** concentration in µg/dL. To express the result in µg/dL calculate with this formula:

$(As/Ast) \times 31.47 = \mu\text{mol/l}$ of Copper.

URINE	R/B	ST _{Urine}	S
Working reagent (R1+R2)	1 (0.5+0.5) mL	1 (0.5+0.5) mL	1 (0.5+0.5) mL
R4 Blank urine	300 µL	---	---
R5 – CAL urine	---	300 µL	---
Sample	---	---	300 µL

Mix carefully. Read the absorbance of Standard (Ast) and of the sample (As) against blank reagent after 10 minutes at 37°C. The colour is stable 30 minutes at room temperature.

CALCULATION URINE

$(As / Ast) \times 10 * = \mu\text{g/dL}$ of Copper.

* is the concentration of **R5- CAL** expressed in µg/dL. To express the result in µg/24 h calculate with this formula:

Copper excreted during 24-hour in µg/24 h =

Copper in µg/dL x 10 x volume of Urine voided in 24-h (in Liters).

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analysers could be completely different by what has been developed as manual determination.

XII. REFERENCE VALUES (see References 1-5)

Normal Values Copper:

Men	70 - 140 µg/dL	(11.0 – 22.0 µmol/L)
Women	80 - 155 µg/dL	(12.6 – 24.4 µmol/L)
Babies (birth-6mths)	20 - 70 µg/dL	(3.15 – 11.0 µmol/L)
Urine	≤ 60 µg/24h	

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

ANALYTICAL PERFORMANCES (Validate on MINDRAY BS300)

The performances of the Reagent **RAME** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 600 µg/dL if it is used with serum/plasma method. For Copper concentrations higher than 600 µg/dL, it is recommended to dilute the sample 1:2 with saline solution, test again and multiply the result x 2.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 8 (1.7) µg/dL in the case of serum / plasma (urine).

Interferences: see References 2.

Interference test criterion: recovery ± 10% of initial value. No interference found on samples with:

- total bilirubin up to 10 mg/dL;
- haemoglobin up to 300 mg/dL,
- lipemia [Intralipid ®] up to 3000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples

The results obtained are following:

Sample	Mean (U/L) ± 2s	CV %
Human 1 urine	20.7 ± 1.7	4.0
Human 2 serum	243.7 ± 3.2	0.7

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (U/L) ± 2s	CV %
Human 1 urine	20.6 ± 2.0	4.4
Human 2 serum	243.8 ± 3.0	0.6

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation $y = 0.9998x + 1.194$
 Correlation coefficient $r = 0.9999$ $n = 20$

XIV. REFERENCES (see Italian version)

 IVD	In Vitro Diagnostic Medical Device	 Temperature limitation	 LOT	Batch code (XXX)	 Manufacturer	 Keep dry	 Non-sterile	Non-sterile
 Consult Instructions for use	Use by (year/month)	 REF	Catalogue number	 Do not reuse	 Fragile, handle with care	 Keep away from heat		

CONTENT

 **NB12000**

R1 - BUFFER	1 x 50 mL
R2 - DiBr-PAESA	1 x 50 mL
R3 - CAL STD 200 serum/plasma	1 x 5 mL
R4 - Blank urine	1 x 15 mL
R5 - CAL STD 10 urine	1 x 5 mL
Instruction for use	1 item

EDMA (EDMS) CODE 11.02.01.06