ADA (Adenosine Deaminase)

IVD

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ADA (Adenosine Deaminase)

IVD

ADENOSINE DEAMINASE QUANTITATIVE UV ASSAY KINETIC METHOD on SERUM, PLEURIC LIQUID and CSF (Cerebrospinal fluid) with MINDRAY BS-300 Analyzer. Traceable to ERM (BCR - 647)

500 tests (5 x 20 mL)

REF NCADA016

INTENDED USE

Adenosine deaminase (ADA) is one of the key enzyme of purine metabolism. It is needed for the breakdown of adenosine from food and for the turnover of nucleic acids in tissues.

Present in virtually all mammalian cells, its primary function in humans is the development and maintenance of the immune system. However, the full physiological role of ADA is not yet completely understood.

The enzyme has been found in bacteria, plants, invertebrates, vertebrates, and mammals, with high conservation of amino acid sequence. The high degree of amino acid sequence conservation suggests the crucial nature of ADA in the purine salvage pathway. However, ADA association has also been observed with epithelial cell differentiation, neurotransmission and gestation maintenance; an increase of ADA activity in pleural fluid is useful in the assessment of tuberculous pleural effusion.

Some mutations in the gene for adenosine deaminase cause it to be not expressed. The resulting deficiency is one cause of lethal severe combined immunodeficiency (SCID).

Deficient levels of ADA improve deoxyadenosine storage, which become shortly dATP, a very toxic inhibitor of the proliferation of leukocytes, killing DNA synthesis; deficient levels of ADA has also been associated with pulmonary inflammation, thymic cell death and defective T-cell receptor signalling.

Elevated levels of ADA has also been associated with AIDS.

PRINCIPLE

The Adenosine Deaminase (ADA) catalyzes the hydrolysis of the substrate Adenosine; the Ammonia obtained from this reaction, reacts with α -ketoglutarate and NADPH by means of Glutamate dehydrogenase (GIDH) giving NADP.

The decrease of absorbance of NADPH, for oxidation to NADP, is proportional to the activity of the ADA in the sample.

Using the standard contained in the kit, traceable to the BCR - 647 International Standard, it is possible to prepare a Calibration Curve to refer.

Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample.

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of ADA higher than 150 U/L, dilute the sample 1:4 with saline solution, repeat the determination and multiply the result by 4.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents).

However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed

the operating procedures for the manipulation of this product.

Material safety data sheet should be supplied on request.

ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of ADA (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) Use sample fresh, as soon as available.

F) **DO NOT USE** haemolyzed samples, because blood cells have very high concentration of ADA.

G) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

H) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

I) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

REAGENTS

Components of the kit:

R1 - BUFFER		REF NCADA016
Good Buffer	>20 mmol/L	5 x 20 mL
R2 - GIDH		5 x 20 mL
Adenosine	>0.1 g/L	
GIDH	>1000 U/L	
NADPH	>0.25 mmol/L	
ketoglutarate	>0,1 mmol/L	
Excipients and stabilizers		
R3 - CAL		1 x 1 mL
Human stabilized serum		
NaN ₃	< 0.1%	

STABILITY: the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8 °C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

ADA (Adenosine Deaminase)

IVD

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- ADA CALIBRATOR for CSF Iyo
- ADA CONTROL Iyo
- ADA CONTROL for CSF Iyo

REF NCADACSF
REF NCADACO1
REF NCADACO2

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - GIDH** with 20 mL of **R1 - BUFFER** (practically all the bottle). Mix gently and let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

STABILITY OF THE WORKING REAGENT

The **WORKING REAGENT** is stable 3 days at 2-8°C.

PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R3 - CAL** with 1 mL of distilled water; don't use saline or PBS, only distilled water!

Mix kindly till complete dissolution. Let the reagent reach the room temperature before use. Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

CAUTION

The Calibrator is from human origin and has to be used as a potential transfer of infective pathologies. This product has been tested and found to be negative for HIV, HCV and HBsAg antibodies by an approved method. Because no test method can offer complete assurance that all infectious agents are absent, it is recommended that this product and all the samples be handled as though capable of transmitting infectious disease.

STABILITY OF CALIBRATOR

The calibrator is stable:

- 5 days is stored at 2-8°C;
- till 3 months at -20°C, fractionated in small volumes. FREEZE only ONE TIME. DO NOT REPEAT FREEZING.

MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Saline solution, distilled water, Calibrators and Controls.

Spectrophotometer or automatic analyzer for Clinical Chemistry.

SAMPLES

• Not haemolysed and not lipemic fresh serum.

• Not haemolysed fresh pleuric liquid.

• Not haemolysed fresh cerebrospinal fluid (CSF).

Samples collection in compliance with CLSI (NCCLS)

(see References 3).

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

ANALYTICAL PROCEDURE ON MANUAL

SPECTROPHOTOMETER

- Wavelength: 340 nm (334-365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 3 minutes
- Sample/reagent: 1/20

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, ST: Calibrator, S: Sample

	S	ST
WORKING REAGENT	1000 µL	1000 µL
Sample	50 µL	----
Calibrator	----	50 µL

Mix carefully and incubate for 4 minutes at 37°C.

Measure the initial absorbances for Calibrator (Ast) and Sample (As).

Repeat the readings for at least 3 times at intervals of one minute. Determine the average of the readings of absorbance/minute for Calibrator (ΔAst/min) and Sample (ΔAs/min).

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analyzers could be completely different by what has been developed as manual determination.

REFERENCE VALUES (see References 4, 5, 6)

Normal Values ADA: 6.8 – 18.2 U/L.

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

CALCULATION

Insert the means found in the following formula:

Adenosine Deaminase (U/L 37°C) =

$$\text{values of Calibrator} \times \frac{\Delta\text{As/min. Sample}}{\Delta\text{Ast/min. Calibrator}}$$

ANALYTICAL PROCEDURE WITH CSF on MANUAL SPECTROPHOTOMETER

PREPARATION OF THE WORKING REAGENT

Is the same of the kit, except for the Calibrator and Control. To calibrate with CSF has to be used a particular ADA Calibrator for CSF at low concentration (**REF** NCADACSF).

To grant the correct test performances we suggest to use ADA Control for CSF (**REF** NCADACO2) (see the relative insert).

- Wavelength: 340 nm (334-365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 5 minutes
- Sample/reagent: 1/4

Let reagents reach the working temperature before use.

ADA (Adenosine Deaminase)

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, ST: Calibrator, S: Sample

	S	ST
WORKING REAGENT	1000 μ L	1000 μ L
Sample	250 μ L	----
Calibrator	----	250 μ L

Mix carefully and incubate for 4 minutes at 37°C.

Measure the initial absorbance for Calibrator (Ast1) and Sample (As1).

Repeat the readings 5 minutes later for Calibrator (Ast2) and Sample (As2).

Calculate the differences of absorbance/5 minutes for

Calibrator Δ Ast/ 5min = Ast1 - Ast2 and

Sample Δ As/ 5min = As1 - As2

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analyzers could be completely different by what has been developed as manual determination.

REFERENCE VALUES (see References 1)

Not available Reference Values for CSF in literature.

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

CALCULATION

Insert the values found in the following formula:

Adenosine Deaminase (U/L 37°C) =

$$\text{values of ADA Calibrator for CSF} \times \frac{\Delta\text{As/ 5min. Sample}}{\Delta\text{Ast/ 5min. Calibrator}}$$

ADDITIONAL NOTE with CSF

1. For Calibration in automation we suggest to apply as follows:

CAL 1 ---- Distilled water (conc. 0.0)

CAL 2 --- Calibrator for CSF (conc. mentioned in Cod. ADACSF)

2. We suggest don't dilute further the Calibrator for CSF.

If absolutely required, REM to dilute the Calibrator for CSF ONLY with distilled water.

USE IMMEDIATELY: stability at 2-8°C no more than half an hour.

ANALYTICAL PERFORMANCES

(validate on MINDRAY BS300)

The performances of the Reagent ADA (Adenosine Deaminase) have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 150 U/L.

However, for ADA concentrations higher than 150 U/L, it is recommended to dilute the sample 1:4 with saline solution, test again and multiply the result x 4.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 1.9 U/L.

Interferences: see References 2.

Interference test criterion: recovery \pm 10% of initial value. No interference found on samples with:

- total bilirubin up to 10 mg/dL;
- haemoglobin up to 300 mg/dL,
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	42.4 \pm 1.8	2.2
Human 2	8.5 \pm 1.4	8.5

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (U/L) \pm 2s	CV %
Human 1	42.9 \pm 2.1	2.5
Human 2	8.9 \pm 1.2	6.6

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation $y = 1.0013x + 0.163$

Correlation coefficient $r = 0.9993$ $n = 20$

REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. CLSI(NCCLS) GP44-A4/H18-A4: Proc. for the Handling and Processing of Blood Specimens for Common Lab. Tests.
4. Martinek R.C., Clin. Chem. 9, 620 (1963).
5. Delia S. et al., Clin. Chem. 33, 1675 (1987).
6. Internal study Maruho Co cod. 5502543.

EDMA (EDMS) CODE 11 01 01 90 00