## CITRIC ACID

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### BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. Tietz Textbook of Clin. Chem. and Molecular Diagnostics, W.B. Saunders Co., Philadelphia (2012), (119; 154; 581).

CODICE RAMO CND W01010299

## CITRIC ACID

CITRIC ACID QUANTITATIVE UV ASSAY on SPERMATIC LIQUID and URINE  
with MINDRAY BS-300 Analyzer

400 tests (6 x 20 mL)

**REF** NCI8822

### INTENDED USE

Citric acid is a weak organic acid and in biochemistry its relative salt, citrate, is an important intermediate in the Krebs cycle.

Pyruvate dehydrogenase complex forms acetyl CoA, from pyruvate and five cofactors then citrate synthase catalyzes the condensation of oxaloacetate with acetyl CoA to form citrate. High citrate conc. can inhibit phosphofructokinase, catalyst of one of the rate-limiting steps of glycolysis. Into the cytoplasm, citrate can also be broken down into acetyl-CoA for fatty acid synthesis. It modulates positively acetyl-CoA carboxylase, enzyme which converts acetyl-CoA into malonyl-CoA (the first step in fatty acid synthesis). It is able to chelate calcium ions, can therefore inhibits coagulation.

The test in urine is used to diagnose renal tubular acidosis and a tendency to form calcium kidney stones. Chronic kidney failure, diabetes, hypoparathyroidism, excessive muscle activity, ACE inhibitors may decrease urine levels.

May increase it an high carbohydrate diet, estrogen therapy, increase of vitamin D.

Citric acid may also be determining citric acid in the sperm, in the ejaculate. In man it is formed in the prostate.

Its concentration in the sperm, as for fructose, depends on the testosterone production of the Leydig's cells; so the citric acid determination appear useful as a functional test for testosterone level in the body and a strong instrument in fertility investigations. Low citric acid values have been found in acute or chronic prostatitis.

### PRINCIPLE

The Citric acid is changed in oxalacetate and acetate by CL (citrate lyase). This reaction is helped by a secondary one, who transforms the oxalacetate originated before (and its decarboxilated product, pyruvate) in the presence of LDH (lactate dehydrogenase), MDH (malate dehydrogenase) and NADH.

The intensity of the UV-colour at this wavelength is proportional to the concentration of Citric acid in the tested sample.

Using the standard contained in the kit **REF** TOG111 it is possible to prepare a Calibration Curve to refer. Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample

### PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of Citric Acid higher than 400 mg/L (2.08 mmol/L), dilute the sample 1:4 with saline solution, repeat the determination and multiply the result by 4.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification,

labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

### ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Citric Acid (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

G) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

### REAGENTS

Components of the kit:

<b>R1 - BUFFER</b>		<b>REF</b> NCI8822
Good buffer	> 10 mmol/L	2 x 70 mL
<b>R2 - MDH, NADH</b>		6 x 20 mL
LDH	> 500 U/L	
NADH	> 0.1 mmol/L	
MDH	> 200 U/L	
<b>R3 - CL</b>		6 x 2 mL
CL	> 300 U/L	
<b>R4 - DILUENT</b>		1 x 12 mL
Good buffer	> 10 mmol/L	
Aminoacids		
Stabilizers		
NaN <sub>3</sub>	< 0.1%	

**STABILITY:** the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8 °C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

## CITRIC ACID

### AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- CI+OX CALIBRATOR liquid
- SUBSTRATE ELEVATED CONTROL Iyo
- SUBSTRATE LOW CONTROL Iyo

REF CTOG111  
REF OG3005  
REF OG3006

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

### PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - MDH, NADH** with 20 mL of **R1 - BUFFER**, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

### STABILITY OF THE WORKING REAGENT

The **WORKING STARTER** is stable:

- 12 days at 2-8°C.;
- 5 weeks at -20°C, when frozen. After reconstitution FREEZE only ONE TIME, DO NOT REPEAT FREEZING.

### PREPARATION OF THE WORKING STARTER

**Rem that this product is temperature sensitive (>30°C).**

Dissolve a vial of **R3 - CL** with 2 mL of **R4 - DILUENT**, mix gently avoid foaming.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

### STABILITY OF THE WORKING STARTER

The **WORKING STARTER** is stable 12 days at 2-8°C.

### MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Saline solution.

Spectrophotometer or automatic analyzer for Clinical Chemistry.

### SAMPLES

- Urine.
- Spermatic liquid.

Samples collection in compliance with CLSI (NCCLS) (see References 3).

The sample can be stored at 2-8°C, up to 6 days.

(see References 1).

### SAMPLES PRETREATMENT

- Use limpid samples; If necessary, filter the samples.
- Use samples at neutral pH; neutralize the Acidic Samples at a neutral pH (max pH ≤ 8.0).
- Proteic samples may be required deproteinization as follows:
  - add perchloric acid (0.33N) 2 mL  
to the sample 0.2 mL
  - mix and centrifuge 5 minutes x 5000 rpm
  - keep 1 mL of the supernatant
  - neutralize with 1 mL sodium hydroxide (NaOH 0.3 N)

- use directly this solution as sample.

- REM that the Dilution Factor will be  $F = 22$ ; keep in mind for the Calculation

- if it will be need further dilutions, please return to PRECAUTIONS FOR USE note 4).

### WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

### ANALYTICAL PROCEDURE ON MANUAL SPECTROPHOTOMETER

- Wavelength: 340 nm (334 - 365 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 8 - 13 minutes
- Sample/reagent: 1/40/4

**Let reagents reach the working temperature before use.**

Pipette in a test tube or cuvette so labelled:

R/B: Reagent Blank, S: Sample

	R/B	S
<b>WORKING REAGENT</b>	1000 $\mu$ L	1000 $\mu$ L
Distilled water	25 $\mu$ L	----
Sample	----	25 $\mu$ L

Mix and incubate for about 3 minutes at 37°C. Measure the absorbance AS1 and AR/B1. Then add:

<b>WORKING STARTER</b>	100 $\mu$ L	100 $\mu$ L

Mix carefully, incubate at 37°C and wait the end of the reaction (5-10 minutes). Read AR/B2, AST2 and AS2.

Calculate for the Reagent/Blank  $AR/B = (AR/B1 - AR/B2)$ ;

Calculate for the Calibrator  $AST = (AST1 - AST2)$ ;

Calculate for the Sample  $AS = (AS1 - AS2)$ .

### ATTENTION!

**The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analysers could be completely different by what has been developed as manual determination.**

### REFERENCE VALUES (see References 1)

Normal Values Citric Acid:

Adults MALES 115 - 922 mg/24 h (0,6 - 4,8 mmol/24 h)

Adults FEMALES 250 - 1153 mg/24h (1,3 - 6,0 mmol/24 h)

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

### CALCULATION

Use this general formula to calculate the concentration:

$(AS - AR/B)$

----- x conc. Calibrator (in mg/L) = mg/L of Citric Acid

$(AST - AR/B)$

### ANALYTICAL PERFORMANCES

(validate on MINDRAY BS300)

The performances of the Reagent **CITRIC ACID** have been tested

## CITRIC ACID

**IVD**

with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

**Method Limitations:** are not know limitations.

**Method Linearity:** the test is linear up to 400 mg/L (2.08 mmol/L). For concentrations  $\geq 400$  mg/L (2.08 mmol/L), it is recommended to dilute the sample 1:4 with saline solution, test again and multiply the result x 4.

**Method Sensitivity (LoD):** the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 8 mg/L (0.042 mmol/L).

**Interferences:** see References 2.

Interference test criterion: recovery  $\pm 10\%$  of initial value. No interference found on samples with:

- total bilirubin up to 40 mg/dL;
- haemoglobin up to 600 mg/dL,
- lipemia [Intralipid®] up to 4000 mg/dL;
- ascorbic acid up to 50 mg/dL.

**Within-run Precision:** determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (mg/L) $\pm 2s$	CV %
--------	----------------------	------

Human 1	212 $\pm$ 8	1.71
Human 2	391 $\pm$ 4	0.55

**Run-to-run Precision:** determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (mg/L) $\pm 2s$	CV %
Human 1	209 $\pm$ 8	1.69
Human 2	398 $\pm$ 6	0.80

**Accuracy:** a group of 20 urine has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation  $y = 0.9894x + 0.005$   
Correlation coefficient  $r = 0.9988$   $n = 20$

### REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. Tietz Textbook of Clin. Chem. and Molecular Diagnostics, W.B. Saunders Co., Philadelphia (2012), (119; 154; 581).