



## INSTRUCTIONS FOR USE

**STREP GROUPING RAPID LATEX TEST KIT**

RAPID LATEX AGGLUTINATION SLIDE TEST FOR GROUPING STREPTOCOCCI OF LANCEFIELD GROUPS A,B,C,D,F AND G FROM CULTURE PLATES

**1 – CLINICAL SIGNIFICANCE AND INTENDED USE**

Most strains of Streptococci which have been isolated from human infections possess group specific antigens. Identification of the organisms includes extraction and characterisation of these antigens from organism grown in culture. The use of specific antisera combined with an enzyme extraction method provide a speedy and reliable typing system for Streptococci.

The streptococcal grouping system provides an enzyme reagent for rapid extraction of the carbohydrate antigens and a series of latex agglutination reagents, specific for groups A, B, C, D, F and G, for rapid detection and identification of the extracted antigens.

Strep Grouping Rapid Latex Test Kit is a manual rapid latex agglutination slide test for grouping Streptococci of Lancefield groups A,B,C,D,F and G from culture plates. The kit is intended for professional use only.

For *in Vitro* diagnostic use only**2 - PRINCIPLE OF THE METHOD**

Latex particles in the Strep Grouping Rapid Latex Test Kit are individually sensitised with rabbit antibodies specific to one of the Streptococcal carbohydrate antigens of groups A, B, C, D, F or G. Streptococcal colonies from culture plates are incubated in an enzyme solution to extract the antigen. The extract / antigen preparation is tested on a reaction card against six suspensions of antibody coated latex particles, each specific for one of the groups A, B, C, D, F and G. In the presence of homologous antigen, particles in one of the suspensions will aggregate to give visible agglutination in contrast to the other suspensions, which will remain un-agglutinated.

**3 - MATERIALS PROVIDED – PACKAGING**

| Product   | Type                     | REF                  | Pack  |
|---|--------------------------|----------------------|---|
| STREP GROUPING RAPID LATEX TEST KIT<br>CND: W0104080302<br>EDMA: 14.02.03.02;<br>RDM: 1555349/R | Latex agglutination test | 271070<br>(50 tests) | <b>REAGENT TEST GR A:</b> 2.5mL – 1 glass bottle contain rabbit Strep Group A antibody-sensitised latex particles in buffer with stabiliser and sodium azide 0.098% as preservative. White cap<br><b>REAGENT TEST GR B:</b> 2.5mL – 1 glass bottle contain rabbit Strep Group B antibody-sensitised latex particles in buffer with stabiliser and sodium azide 0.098% as preservative. White cap<br><b>REAGENT TEST GR C:</b> 2.5mL - 1 glass bottle contain rabbit Strep Group C antibody-sensitised latex particles in buffer with stabiliser and sodium azide 0.098% as preservative. White cap<br><b>REAGENT TEST GR D:</b> 2.5mL - 1 glass bottle contain rabbit Strep Group D antibody-sensitised latex particles in buffer with stabiliser and sodium azide 0.098% as preservative. White cap<br><b>REAGENT TEST GR F:</b> 2.5mL – 1 glass bottle contain rabbit Strep Group F antibody-sensitised latex particles in buffer with stabiliser and sodium azide 0.098% as preservative. White cap<br><b>REAGENT TEST GR G:</b> 2.5mL – 1 glass bottle contain rabbit Strep Group G antibody-sensitised latex particles in buffer with stabiliser and sodium azide 0.098% as preservative. White cap<br><b>CONTROL +:</b> 1.0mL – 1 glass bottle of Positive control: contains inactivated polyvalent antigenic extracts to groups A,B,C,D,F and G preserved with 0.098% sodium azide. Red cap<br><b>ENZ:</b> 2 x 10mL – 2 glass bottles with lyophilised extraction enzyme<br>Slide, 6 test areas: plastic waterproof sheets for reaction (50 items)<br>Sticks (1x25): plastic sticks for mixing (12 items)<br>Secondary packaging: cardboard box. |

**4 - MATERIALS REQUIRED BUT NOT PROVIDED**

Bacteriological loops. Glass or plastic test tubes. Pipette to dispense 0.4 mL volume. Water bath set at 37°C. Sample droppers or Pasteur pipettes. Timer or clock.

**5 - PRECAUTIONS AND WARNINGS**

- STREP GROUPING RAPID LATEX TEST KIT is for *in vitro* diagnostic, for professional use only; it is to be used by adequately trained and qualified laboratory personnel.
- The sensitivity of the test may be reduced at low temperatures. Allow the reagents and samples to reach room temperature (15-30°C/59-86°F) before use.
- The slides are made of plastic, washed with distilled water. If the test area of the slide becomes water resistant, clean it with alcohol.
- Do not use after expiration date or if the packaging is damaged. The quality of the reagent cannot be guaranteed beyond their shelf-life date or if the reagents are stored under inappropriate conditions.
- Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Appropriate precautions should be taken when handling or disposing of potential pathogens. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30minutes. Liquid waste containing acid must be neutralised before treatment.
- All operations linked to the use of the test must be performed in accordance with Good Laboratory Practices and according to the Instruction for use of the kit.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The Certificates of Analysis and the Safety Data Sheet of the product are available on the website [www.masciabrunelli.it](http://www.masciabrunelli.it).
- The information provided in this document has been defined to the best of our knowledge and ability and represents a guideline for the proper use of the product but without obligation or liability. In all cases existing local laws, regulations and standard procedures must be observed for the examination of samples collected from human and animal organic districts, for environmental samples and for products intended for human or animal consumption. Our information does not relieve our customers from their responsibility for checking the suitability of our product for the intended purpose.
- Notify Mascia Brunelli Spa and the Relevant Authorities of any serious incidents occurring in connection with the *in vitro* diagnostic device. [complaint@masciabrunelli.it](mailto:complaint@masciabrunelli.it)

**6 - STORAGE CONDITIONS AND SHELF LIFE**

All the kit components will remain stable until the expiration date printed on the label, *when stored tightly closed at 2-8°C and contaminations are prevented during their use*. Do not freeze: *frozen reagents could change the functionality of the test*. Extraction Enzyme is stable for 3 months after reconstitution if stored at 2-8°C. To prolong the life of the enzyme, it may be dispensed into suitable test tubes in 0.4mL volumes and stored frozen, at -20°C or below when it will be stable for 6 months. **Enzyme should not be frozen and thawed more than once.**





Deterioration of reagents should be suspected if:

- Clumping of any of the latex reagents is evident and cannot be removed by shaking vigorously for a few seconds.
- The Positive Control or Extraction Enzyme becomes cloudy or forms a sediment.
- The Positive Control fails to cause agglutination of one or more latex reagents within the recommended reaction time.
- Un-inoculated Extraction Enzyme causes agglutination of any of the latex reagents.

Reagents showing signs of deterioration should not be used.

### 7 – SPECIMENS COLLECTION

The normal media used for culture preparations include blood agar base, in such cases note colonial characteristics, haemolysis, and cell morphology prior to testing. Ensure the organisms to be tested are Gram-positive and catalase negative. Any blood agar plate culture with 2-6 separate colonies may be used, they should have been inoculated from a pure culture of the organism. If a conclusive result of cultures that appear to contain streptococci is not obtained, further subculture of suspect colonies is recommended. Organisms of groups A, B, C, D, F and G are normally beta haemolytic. Any alpha or non-haemolytic organisms showing positive results should be confirmed by further biochemical tests. (Some group B and D strains can be either alpha or non-haemolytic).

### 8 - TEST PROCEDURE

#### Quality Control

The Positive Control should be tested regularly to ensure that the reagents are functioning correctly.

The control is supplied ready for use and should be tested in place of the culture extract in the test procedure. The Positive Control should give agglutination with all the test Latex Reagents. Failure of the Positive Control to give an agglutination pattern may be evidence of latex reagent deterioration. If a negative control is desired, un-inoculated Extraction Enzyme should be tested in place of the culture extract in the test procedure. **Reactions containing traces of indistinct granulation/agglutination may be observed; these should be ignored and considered negative.**

#### Test Procedure

Proceed as follows for each organism to be grouped.

1. **Allow the Latex reagents and positive control to reach room temperature.**
2. Just prior to use, reconstitute a bottle of enzyme by adding 10mL distilled water. Mix gently to ensure complete reconstitution. Dispense 0.4mL **Extraction Enzyme** into a test tube.
3. Pick Streptococcal colonies from the surface of the agar using a bacteriological loop and emulsify them thoroughly in the Extraction Enzyme. To obtain best results, pick **2-4 colonies (2-3 mm diameter)** or their equivalent for extraction. Excessive inoculation of extraction enzyme may cause non-specific agglutination. For minute-colony strains, a sweep of growth will be necessary (**if a broth culture is to be grouped, pipette 0.1 mL of an overnight culture into 0.4 mL extraction enzyme**).
4. Incubate the tube for **10 minutes in a 37°C water bath**. Shake the tube after the first 5 minutes incubation and shake vigorously prior to testing to obtain even suspension of antigen.
5. Vigorously shake Latex reagents for a few seconds to obtain even suspension. Dispense one drop of each **Latex reagent** separately into six circles on a reaction card.
6. Transfer one drop of **well mixed extract** (or Positive Control) into the six separate circles next to the drop of latex reagent.
7. Mix the contents of each circle using separate mixing sticks and spread the liquid to cover the area of the circle. Do not use the same mixing stick for more than one circle.
8. Slowly and gently, rock and rotate the reaction card to mix the reagents for a maximum of **one minute**.
9. Inspect the card for agglutination. If present, agglutination should be clearly visible with the naked eye.

### 9 – READING, INTERPRETATION

**POSITIVE RESULT:** Indicated by rapid strong aggregation of the latex particles with one group reagent (Figure 1). Subsequent reactions in other circles with the same extract should be disregarded. Only strong agglutination is significant. This will normally occur within a few seconds of mixing; however, the time is dependent on the extract strength.

**NEGATIVE RESULT:** Indicated by a milky appearance, without any significant aggregation of the latex particles (Figure 2). Traces of indistinct aggregation should be ignored and considered negative.

**INCONCLUSIVE:** With weaker extracts agglutination may take longer than 1 minute to appear and give smaller clumps. If this occurs tests should be repeated with a fresh subculture. If the same result is observed after retesting, alternative biochemical methods should be conducted to identify the isolate.

**NON-SPECIFIC RESULT:** Occasionally, strains of streptococci may give weak reactions with more than one group. If this occurs tests should be repeated. If agglutination occurs in all groups, either the extraction enzyme has been over-inoculated in which case repeat the test using a lighter inoculum, or a mixed culture was tested, in which case subculture and retest. Boiling the remaining extract for two or three minutes, cooling and retesting may lead to clearer results.

Figure 1

Figure 2



**Figure 1:** a positive result is indicated by visible aggregation of the latex particles.

**Figure 2:** a negative result is indicated by a milky appearance without any visible aggregation of the latex particles.

Sensibilità = 92% (607/662)

Specificità = 100% (24/24)

### 10 – NOTES

#### Colonies associated with beta-haemolysis:

1. Agglutination of a single latex reagent indicates the group identity of the strain. Complimentary tests should be considered to confirm the results, in particular:
  - for group D strains, biochemical tests to differentiate *Enterococcus* species from group D *streptococcus* species (the former has relatively high antibiotic resistance).
  - for group A, C or G strains with minute colony morphology, biochemical tests to confirm *S.millieri* / *S.anginosus* identification.
2. Agglutination of more than one latex reagent indicates the possibility of mixed growth of organisms from different groups or the presence of a strain with more than one group (for example some group D streptococci which also possess group G antigen).





Further procedures to be considered:

- subculture to obtain pure isolates for retesting.
- for strains with group D and group G antigen, biochemical tests to differentiate Enterococcus species from group D streptococcus species (Enterococcus strains with both these antigens may be more antibiotic resistant than those with only group D antigen).

3. No significant agglutination in any of the latex reagents indicates either that no group A, B, C, D, F or G streptococci were present in the test sample or that they were present in numbers below the threshold of sensitivity of the test.

Further procedures to be considered:

- retest using a higher inoculum, particularly if group D or group F streptococci are suspected.
- beta-haemolytic streptococci which do not group may be

**Colonies not associated with beta-haemolysis:**

Agglutination of a single latex reagent showing a result of group B or group D gives a reliable identification of the strain. If the result is group A, C, F or G it may not be relevant to the identification of the strain and other identification methods are necessary.

Further procedures to be considered:

- if the result is group D, biochemical differentiation between Enterococci and group D streptococci (see above).

Any other combination of results should be interpreted using the information provided above.

**11 – CHARACTERISTICS NOTES**

Strep Grouping Rapid Latex Test Kit has been evaluated in comparison with a leading commercial latex kit as a reference for grouping Streptococci, using clinical samples at a number of independent sites. Overall Results are shown in Table.

Sensitivity: 607/662 = 92%  
Specificity: 24/24 = 100%

|                        |     | Strep Grouping Rapid Latex Test Kit |     | Total |
|------------------------|-----|-------------------------------------|-----|-------|
|                        |     | +ve                                 | -ve |       |
| Leading commercial kit | +ve | 607                                 | 55  | 662   |
|                        | -ve | 0                                   | 24  | 24    |
| Total                  |     | 607                                 | 79  | 686   |

**REPRODUCIBILITY**

**Intra Batch reproducibility** was evaluated by testing sensitive ty of one batch of each of the test latexes on ten separate occasions with three different operators against serial dilutions of reference antigens. End point titres varied by a maximum of one doubling dilution from assay to assay.

**Inter Batch Reproducibility** was examined by testing sensitivity and specificity of 10 batches of product against serial dilutions of reference antigens. Between the batches variation in titres was a maximum of one.

**12 – LIMITATIONS OF THE METHOD**

- As with all diagnostic tests, a final diagnosis cannot rely on the outcome of a single test but the results should be interpreted in the context of all available clinical and laboratory information.
- Accurate results depend on testing an appropriate amount of growth.
- This is not usually a problem, however some strains of streptococci belonging to group D possess lower or negligible quantities of group antigen and some strains of group F may be difficult to remove from the surface of agar plates. Antigen production in group D strains may be improved by culturing them on agar supplemented with 0.5 to 1.0% glucose. This supplement does obscure demonstration of haemolysis but it may be considered in situations where it is important to resolve identification.
- Growth of minute-colony strains may be improved by culture in a carbon dioxide enriched atmosphere.
- Streptococci from groups Q, R and S may also possess detectable levels of group D antigen.
- Antigens common to the streptococcal group antigens have been described in a number of unrelated species. For example false positive reactions can occur with *Escherichia*, *Klebsiella* or *Pseudomonas*. These are normally easily differentiated by cultural characteristics and cause no confusion in streptococcal identification.
- The components of this I.v.D. were always tested together without compatibility with components from other manufacturers. While not excluding the possibility that these components can be used with components of the same formulation but produced by other companies, there is no experimental evidence of such compatibility.

**13 – REFERENCES**

1. Birch, B.R. et al (1984). Antibiotic susceptibility and biochemical properties of *Streptococcus faecalis* strains reacting with both D and G antisera. *J. Clin Pathol* 37, 1289-1292.
2. Elliot, S.D. and Tai, J.Y. (1978). The type-specific polysaccharides of *Streptococcus suis*. *J Exp Med* 148, 1699-1704.
3. Facklam, R.R. (1977). Physiological differentiation of viridans streptococci. *J Clin Microbiol* 5, 184-201.
4. Facklam, R.R. (1985). Serologic identification of streptococci : how useful is serologic grouping? *Clin Microbiol Newsletter* 7, 91-94.
5. Facklam, R.R. and Smith, P.B. (1976). The Gram positive cocci. *Human Pathology* 7, 187-194.
6. Facklam, R.R. and Washington, J.A. (1991). Streptococcus and related catalase-negative Gram-positive cocci. In manual of Clinical Microbiology 5th Ed, Edited by Balows, A et al American Society for Microbiology, Washington D.C. Pages 238-257.
7. Harvey, C.L. and McIlmurray, M.B. (1984). Streptococci with dual antigen specificity for Lancefield groups D and G. *Eur j Clin Microbiol* 3, 526-530.
8. Lancefield, R.C. (1933). A serological differentiation of human and other groups of haemolytic streptococci. *J. Exp med* 57, 571-595.
9. Medrek, T.F. and Barnes, E.M. (1962). The influence of the growth medium on the demonstration of a group D antigen in faecal streptococci. *J Gen microbial* 28, 701-709.
10. Nowlan, S.S. and Deibel, R.H. (1967). Group Q streptococci 1. Ecology, serology, physiology and relationship to established enterococci. *J. Bacteriol* 94, 291-296.
11. Parker, M.T. and Ball, L.C. (1976). Streptococci and Aerococci associated with infection in man. *J. Med Microbiol* 9, 275-302.

**TABLE OF APPLICABLE SYMBOLS**

|  |                                    |  |                        |  |                   |  |              |  |                           |  |                          |
|--|------------------------------------|--|------------------------|--|-------------------|--|--------------|--|---------------------------|--|--------------------------|
|  | In Vitro Diagnostic Medical Device |  | Temperature limitation |  | Batch code (DXXX) |  | Manufacturer |  | Keep dry                  |  | Unique device identifier |
|  | Consult Instructions for use       |  | Use by (year/month)    |  | Catalogue number  |  | Do not reuse |  | Fragile, handle with care |  | Keep away from heat      |

**REVISION HISTORY**

| Version                                 | Description of changes     | Date    |
|---|----------------------------|---------|
| Instructions for Use (IFU) - Revision 3 | Updated layout and content | 2022/09 |

Note: minor typographical, grammatical, and formatting changes are not included in the revision history.

