

**INSTRUCTIONS FOR USE****MULLER KAUFFMANN TETRATHIONATE BROTH BASE**

Dehydrated culture medium



Muller Kauffmann Tetrathionate Broth – from the left: uninoculated tube and tube with growth of *S. Typhimurium*.

1 - INTENDED USE

In vitro diagnostic. Selective liquid medium for the enrichment of *Salmonella* from food and faecal specimens.

**2 - COMPOSITION - TYPICAL FORMULA*
(AFTER RECONSTITUTION WITH 1 L OF WATER)**

Tryptone	7.00 g
Soy Peptone	2.30 g
Sodium Chloride	2.30 g
Calcium Carbonate	25.00 g
Sodium Thiosulphate	40.70 g
Bile Salts	4.75 g

*The formula may be adjusted and/or supplemented to meet the required performances criteria.

3 - PRINCIPLE OF THE METHOD AND EXPLANATION OF THE PROCEDURE

Muller Kauffmann Tetrathionate Broth has been originally described by Muller¹ and later modified by Kauffmann² by the inclusion of ox bile and brilliant green as selective agents to suppress bacteria such as *Proteus* spp.

The medium is included in the reviews ISTISAN 05/27³ and ISTISAN 96/35⁴ for the selective enrichment of *Salmonella* from samples of the food chain and from faeces prior to selective isolation.

Tryptone and soy peptone provide carbon, nitrogen, vitamins and minerals for microbial growth; the selective agents of the medium are bile salts, the added brilliant green and sodium tetrathionate which is formed from the sodium thiosulfate when the iodine / potassium iodide solution is added to the medium; calcium carbonate neutralizes the sulfuric acid that is produced by the reduction of tetrathionate during the growth of salmonellae, keeping the pH at neutral values. The complete medium allows the development of salmonellae and is inhibitory for Gram-positive bacteria and for a large part of Gram-negative bacteria of enteric origin.

4- DIRECTIONS FOR MEDIUM PREPARATION

Suspend 82 g in 1000 mL of cold purified water and bring with frequent agitation. Cool to 42-45°C and, immediately before use, add 19 mL of Jodine Solution (REF 421501) and 9.5 mL of Brilliant Green 0.1% Solution (REF 421505). Mix well and distribute into test sterile tubes (10 mL) or flasks (100 mL).

5 - PHYSICAL CHARACTERISTICS

Dehydrated medium appearance	pale green, fine, homogeneous, free-flowing powder
Solution and prepared plates appearance	pale green with white precipitate
Final pH at 20-25°C	8.0 ± 0.2

6 - MATERIALS PROVIDED - PACKAGING

Product	Type	REF	Pack
Muller Kauffmann Tetrathionate Broth Base	Dehydrated culture medium	4017432	500 g (6,1 L)
		4017434	5 kg (61 L)

7 - MATERIALS REQUIRED BUT NOT PROVIDED

Sterile loops, needles and swabs, incubator and laboratory equipment as required, ancillary culture media and reagents for the identification of the colonies.

8 - SPECIMENS

Muller Kauffmann Tetrathionate Broth may be used for the enrichment of faecal specimens. Collect stool according to standard procedures, with swab preferably with transport medium or in a stool container with or without transport fluid. Good laboratory practices for collection, transport and storage of the clinical specimens should be applied

Food samples: refer to applicable Standards and laws.

9 - TEST PROCEDURE

Allow tubes to come to room temperature. For faeces examination, inoculate test tubes with 1 g of faeces, or 1 mL of faecal suspension obtained suspending 1 g of faeces in 1 mL of saline solution. Rectal swabs received fresh or in transport medium should be rinsed thoroughly in 1 mL of saline.

Incubate the inoculated tubes in aerobic atmosphere at 35-37°C for 18-24 hours.

For milk and dairy products, the following procedure can be used:

Transfer 25 g of sample to 225 mL of Buffered Peptone Water and incubate at 35-37°C for 18-24 hours.





From the pre-enrichment broth transfer 2 aliquots of 10 mL respectively into 100 mL of Muller Kauffmann Tetrathionate Broth and into 100 mL of Selenite Cystine Broth.
Incubate Muller Kauffmann Tetrathionate Broth at 42-43°C for 24 and 48 hours and Selenite Cystine Broth at 35-37°C for 24 and 48 hours. After 24 and 48 hours of incubation subculture on appropriate selective enteric media.

10 - READING AND INTERPRETATION

After incubation, growth of organisms is indicated by turbidity and discolouration. Subculture by streaking a loopful of broth on selective enteric plating media. The plating media should be chosen as a combination of greater and lesser inhibitory selective agars.

11 - USER QUALITY CONTROL

All manufactured lots of the product are released for sale after the Quality Control has been performed to check the compliance with the specifications. However, the end user can perform its own Quality Control in accordance with the local applicable regulations, in compliance with accreditation requirements and the experience of the Laboratory. Here below are listed some test strains useful for the quality control of medium.

CONTROL STRAINS	INCUBATION T°/ t / ATM	EXPECTED RESULTS
S.Typhimurium ATCC 14028	35-37°C / 18-24h / A	good growth after subculture to TSA plate
E.coli ATCC 25922	35-37°C / 18-24h / A	scanty growth after subculture to TSA plate

A: aerobic incubation; ATCC is a trademark of American Type Culture Collection

12 - PERFORMANCES CHARACTERISTICS

Prior to release for sale a representative sample of all lots of dehydrated Muller Kauffmann Tetrathionate Broth Base REF 401743, supplemented with brilliant green and iodine solutions, are tested for productivity and selectivity by comparing the results with a previously approved Reference Batch.

Productivity is tested by dilution to extinction method, by inoculating 1 mL of appropriate decimal dilutions of target organisms in test tubes and incubating at 35-37°C for 18-24 hours and recording the highest dilution showing growth in Reference Batch (G_{RB}) and in Test Batch (G_{TB}) after sub-culture on Tryptic Soy Agar plates. Productivity is tested with the following target strains: S.Typhimurium ATCC 14028, S.Enteritidis ATCC 13076. The productivity index $G_{RB}-G_{TB}$ for each test strain shall be ≤ 1 .

Selectivity is evaluated with dilution to extinction method, by inoculating 1 mL of appropriate decimal dilutions of non-target organisms in duplicate test tubes and incubating at 35-37°C for 18-24 hours and recording the highest dilution showing growth in Reference Batch (G_{RB}) and in Test Batch (G_{TB}) after sub-culture on Tryptic Soy Agar plates. Selectivity is tested with the following non-target strains: E.coli ATCC 25922, E.faecalis ATCC 29212. The selectivity index $G_{RB}-G_{TB}$ for each test strain shall be ≥ 1 .

Productivity and selectivity are tested also together with a mixture of appropriate dilutions of target and non-target strains: S.Typhimurium ATCC 13076+E.coli ATCC 25922+P.aeruginosa ATCC 27853. After incubation of inoculated tubes at 35-37°C for 18-24 hours and subculture on XLD Agar plate, the target strain will show a predominant growth on plated medium.

13 - LIMITATIONS OF THE METHOD

- Muller Kauffmann Tetrathionate Broth is not suitable for growth of S.Typhi, S.Paratyphi, S.Sendai, S.Gallinarum; it is not recommended for examination of typhoid fever.⁵
- After the enrichment in Muller Kauffmann Tetrathionate Broth, even if the microbial colonies on the plates are differentiated on the basis of their morphological and chromatic characteristics, it is recommended that biochemical, immunological, molecular, or mass spectrometry testing be performed on isolates, from pure culture, for complete identification. If relevant, perform antimicrobial susceptibility testing.
- This culture medium is intended as an aid in the diagnosis of infectious diseases; the interpretation of the results must be made considering the patient's clinical history, the origin of the sample and the results of other diagnostic tests.

14 - PRECAUTIONS AND WARNINGS

- This product is a qualitative *in vitro* diagnostic, for professional use only; it is to be used by adequately trained and qualified laboratory personnel, observing approved biohazard precautions and aseptic techniques.
- Dehydrated media must be handled with suitable protection. Before use, consult the Safety Data Sheet.
- Apply Good Manufacturing Practice in the production process of prepared media.
- This culture medium contains raw materials of animal origin. The *ante* and *post mortem* controls of the animals and those during the production and distribution cycle of the raw materials, cannot completely guarantee that this product doesn't contain any transmissible pathogen. Therefore, it is recommended that the culture medium be treated as potentially infectious, and handled observing the usual specific precautions: do not ingest, inhale, or allow to come into contact with skin, eyes, mucous membranes. Download the TSE Statement from the website www.biolifeitaliana.it, describing the measures implemented by Biolife Italiana for the risk reduction linked to infectious animal diseases.
- All laboratory specimens should be considered infectious.
- The laboratory area must be controlled to avoid contaminants such as culture medium or microbial agents.
- Sterilize all biohazard waste before disposal. Dispose the unused medium and the sterilized media inoculated with samples or microbial strains in accordance with current local legislation.
- Do not use the culture medium as active ingredient for pharmaceutical preparations or as production material intended for human and animal consumption.
- The Certificates of Analysis and the Safety Data Sheet of the product are available on the website www.biolifeitaliana.it.
- Notify Biolife Italiana Srl (complaint@biolifeitaliana.it) and the relevant Authorities of any serious incident occurring in connection with the use of the *in vitro* diagnostic.
- The information provided in this document has been defined to the best of our knowledge and ability and represents a guideline for the proper use of the product but without obligation or liability. In all cases existing local laws, regulations and standard procedures must be observed for the examination of samples collected from human and animal organic districts, for environmental samples and for products intended for human or animal consumption. Our information does not relieve our customers from their responsibility for checking the suitability of our product for the intended purpose.





15 - STORAGE CONDITIONS AND SHELF LIFE

Upon receipt, store at +10°C /+30°C away from direct light in a dry place. If properly stored, it may be used up to the expiration date. Do not use beyond this date. Avoid opening the bottle in humid places. After use, the container must be tightly closed. Discard the product if the container and/or the cap are damaged, or if the container is not well closed, or in case of evident deterioration of the powder (colour changes, hardening, large lumps).

The user is responsible for the manufacturing and quality control processes of prepared media and for the validation of the shelf life of the finished products, according to the type (tubes/bottles) and the storage method (temperature and packaging).

16 – REFERENCES

1. Muller, L. (1923) C.R. Soc. Biol. (Paris) **89**, 434-443
2. Kauffmann, F. (1935) Z.f. Hyg. **117**, 26-32
3. Rapporto ISTISAN 05/27. ISSN 1127-3117. Infezioni da Salmonella: diagnostica, epidemiologia e sorveglianza. Raccolta a cura di C.Graziani, P.Galetta, L.Busani, AM Dionisi, E.Filetici, A.Ricci, A.Caprioli, I.Luzzi.
4. Rapporto ISTISAN 96/35. ISSN 1123-3117. Metodi di analisi per il controllo microbiologico degli alimenti. Raccolta a cura di D. De Medici, L. Fenicia, L. Orefice e A.Stacchini.
5. MacFaddin JF. Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria. Baltimore: Williams & Wilkins; 1985.

401743 MULLER KAUFFMANN TETRATHIONATE BROTH BASE

SDS rev 1
Regulation (EU) 2020/878

Contains: BILE SALT

Classification:

Skin sensitization, category 1 H317 May cause an allergic skin reaction.

Labelling

Hazard pictograms:



Signal words: Warning

Hazard statements:

H317 May cause an allergic skin reaction.

Precautionary statements:

- P280** Wear protective gloves.
- P261** Avoid breathing dust / fume / gas / mist / vapours / spray.
- P333+P313** If skin irritation or rash occurs: Get medical advice / attention.
- P362+P364** Take off contaminated clothing and wash it before reuse.

TABLE OF APPLICABLE SYMBOLS

REF or REF Catalogue number	LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	Manufacturer	This side up	Store in a dry place
Temperature limitation	Content sufficient for <n> tests	Consult instructions for Use	Use by	Fragile	Keep away from direct light

REVISION HISTORY

Version	Description of changes	Date
Revision 1	Updated layout and content	2022/03
Revision 2	Removal of obsolete classification	2023/04
Revision 3	Addition of SDS indications	2024/02

Note: minor typographical, grammatical, and formatting changes are not included in the revision history.

