N °541308 BE -0 08/2004 page 1/1

CLOSTRIDIUM DIFFICILE SELECTIVE MEDIUM

Ready to use plates for the isolation of *Clostridium difficile*.

Typical formula (g/l)

Peptocomplex	40.0
Disodium Hydrogen Phosphate	5.0
Potassium Dihydrogen Phosphate	1.0
Magnesium Sulphate	0.1
Sodium Chloride	2.0
Fructose	6.0
Agar	15
D-Cycloserine	500 mg
Cefoxitin	16 mg
Defribinated Horse Blood	50 ml

Final pH 7.4 ± 0.2

Description

Clostridium difficile was first isolated from meconium and infant faeces in 1935 by Hall and O'Toole but only recently has been seen a clear correlation between the cytotoxin produced by Clostridium difficile and the illness characterised by pseudomembranous colitis with inflammatory plaques, serious diarrhoea, dehydration, fever and leucocytosis. The most suitable medium for the selective isolation of Clostridium difficile is that proposed by George et al. in 1978. Biolife's medium is a slightly modified version of the original formula. It contains fructose, cycloserine, cefoxitin and defibrinated horse blood as a substitute for the egg yolk emulsion, used by George et al. Cycloserine and cefoxitin inhibit enterobacteria, faecal streptococci, anaerobic non spore-forming Gram-negative bacilli and Clostridium spp. except for Clostridium difficile.

Technique

Streak the specimen directly onto the agar surface, to obtain well isolated colonies.

Incubate at 37°C for 48 hours in anaerobic conditions.

After 48 hours *Clostridium difficile* colonies are visible with the following characteristics: circular, raised, opaque grey, sometimes with irregular borders, 4-6 mm in diameter.

Storage

Store at 2-8° - When stored as directed the plates remain stable until the expiry date shown on the label. Do not use beyond stated expiry date.

References

- Bartlett, J.G. et al. (1978) N. Eng. J. Med., 298, 531.
- Boriello S.P. et al. (1981) J. Antimicrob. Chemother. 7 Supp. A. 53-62.
- George R.H. et al. (1976) J. Clin. Microbiol. 6, 214-219.
- George, R.H. et al. (1978) Brt. Med. J., i. 695.
- Hall and E. O'Toole (1935) Am. J. Dis. Child. 49, 390.
- Larson, H.E. et al. (1978) Lancet, i. 1063-1066

Packaging

541308 C. Difficile Selective Medium 20 ready to use plates