

**INSTRUCTIONS FOR USE****SABOURAUD DEXTROSE AGAR +CAF+CYCLOHEXIMIDE****Ready-to-use plates**

Trichophyton mentagrophytes
on Sabouraud Dextrose Agar+ CAF+ Cycloheximide

1 - INTENDED USE

In vitro diagnostic device. Selective medium for the isolation of pathogenic fungi, especially dermatophytes, from clinical specimens.

2 - COMPOSITION - TYPICAL FORMULA *

Pancreatic digest of casein	5.00 g
Peptic digest of meat	5.00 g
Glucose	40.00 g
Agar	15.00 g
Chloramphenicol	0.05 g
Cycloheximide	0.50 g
Purified water	1000 mL

*the formula may be adjusted and/or supplemented to meet the required performances criteria.

3 - PRINCIPLE OF THE METHOD AND EXPLANATION OF THE PROCEDURE

By the end of the 1890's, Raymond Jacques Sabouraud had crystallized and organized the scattered observations regarding the role of pathogenic fungi in dermatophytic infections and proposed a medium for their isolation and classification.^{1,2} This formulation was named Sabouraud medium and is the basic routine culture medium used to grow fungi in clinical laboratories.

The components of the basal medium Sabouraud Dextrose Agar conform to the recommendations of the current European Pharmacopoeia³. The addition of chloramphenicol and cycloheximide is a modification designed to increase the selective properties and to improve the isolation of pathogenic fungi, especially dermatophytes, from specimens contaminated with saprophytic fungi and bacteria. Following the initial report of Whiffen, *et al.*,⁴ cycloheximide has been found to be of value in increasing the number of isolations of pathogenic fungi from clinical materials.⁵

Sabouraud Dextrose Agar+CAF+Cycloheximide (SDA CAF-CEX) is particularly useful for the investigation of dermatological specimens for superficial mycosis and the target organisms are dermatophytes and some yeasts.⁶

Pancreatic digest of casein and peptic digest of animal tissue provide nitrogen, carbon and trace elements for microbial growth. The low pH is favourable for the growth of fungi and is slightly inhibitory to contaminating bacteria; glucose, at high concentration is a carbon and energy source. The selective properties are increased by the presence of chloramphenicol, a broad-spectrum antibiotic, which is inhibitory to a wide range of Gram-negative and Gram-positive bacteria and cycloheximide that inhibits the faster-growing saprophytic fungi.⁷

4 - PHYSICAL CHARACTERISTICS

Prepared plates appearance	yellow, limpid
Final pH at 20-25 °C	5.6 ± 0.2

5 - MATERIALS PROVIDED - PACKAGING

Product	Type	REF	Pack
Sabouraud Dextrose Agar+CAF+Cycloheximide	Ready-to-use plates	542008	2 x 10 plates ø 90 mm primary packaging: 2 cellophane sachets secondary packaging: cardboard box

6 - MATERIALS REQUIRED BUT NOT PROVIDED

Sterile loops and swabs, incubator and laboratory equipment as required, ancillary culture media and reagents for the identification of the colonies.

7 - SPECIMENS

Sabouraud Dextrose Agar+CAF+Cycloheximide can be directly inoculated with clinical specimens collected from sites contaminated with saprophytic fungi and bacteria, mainly skin, nail, hair. Consider that cycloheximide may inhibit some opportunistic fungi (see Limitations of the method). Refer to the quoted literature for specimen types, related to specific infections.^{6,8} Sabouraud Dextrose Agar+CAF+Cycloheximide is not suitable for direct inoculation of specimens from normally sterile sites. Collect specimens before antimicrobial therapy where possible. Good laboratory practices for collection, transport and storage of the clinical specimens should be applied; consult appropriate references for further information.^{6,8}

8 - TEST PROCEDURE

Allow plates to come to room temperature. Inoculate the clinical specimen as soon as possible after collection; streak with a loop over the four quadrants of the plate to obtain well isolated colonies. Alternatively, if the material is being cultured directly from a swab, roll the swab over a small area of the surface at the edge; then streak from this inoculated area. For cutaneous samples, press specimen lightly into medium.

The user is responsible for choosing the appropriate incubation time and temperature depending on the processed specimen, the requirements of organisms to be recovered and the local applicable protocols.

For dermatophytes detection incubate at 26-30°C and examine cultures every 4-6 days for a period of up to 21 days.⁶





9 - READING AND INTERPRETATION

After incubation, observe the bacterial growth and record the specific morphological and chromatic characteristics of the colonies and sub-culture to appropriate media for further identification tests.

10 - USER QUALITY CONTROL

All manufactured lots of the product are released for sale after the Quality Control has been performed to check the compliance with the specifications. However, the end user can perform its own Quality Control in accordance with the local applicable regulations, in compliance with accreditation requirements and the experience of the Laboratory. Here below are listed some test strains useful for the quality control.⁹

CONTROL STRAINS	INCUBATION T° / T / ATM	EXPECTED RESULTS
<i>S.cerevisiae</i> ATCC 9763	26-28°C / 72 H / A	inhibited
<i>T.mentagrophytes</i> ATCC 9533	26-28°C / 72 H / A	good growth, white colonies with typical morphology
<i>E.coli</i> ATCC 25922	26-28°C / 72 H / A	inhibited

A: aerobic incubation; ATCC is a trademark of American Type Culture Collection

11 - PERFORMANCES CHARACTERISTICS

Prior to release for sale a representative sample of all lots of ready to use plates of Sabouraud Dextrose Agar+CAF+Cycloheximide is tested for productivity and selectivity.

The productivity characteristics are tested by semi-quantitative ecometric technique with the following target strains: *T.mentagrophytes* ATCC 9533 and *C.albicans* ATCC 10231. After incubation at 26-28°C for up 72 hours, the amount of growth and colonies' characteristics are evaluated and recorded. The target strains show good growth with white colonies and typical morphology.

The selectivity is evaluated with modified Miles-Misra surface drop method by inoculating the plates with suitable decimal dilutions in saline of a 0.5 McFarland suspension of the non-target strains *E.coli* ATCC 25922, *A.brasiliensis* ATCC 16404 and *S.cerevisiae* ATCC 9763. The growth of non-target strains is totally inhibited.

12 - LIMITATIONS OF THE METHOD

- Cycloheximide may inhibit some important opportunistic fungi such as *Fusarium*, *Scopulariopsis*, *Pseudallescheria*, *Trichosporon*, some *Aspergillus* spp., *Talaromyces* (formerly *Penicillium*) *marneffeii*, mucoraceous fungi, some dematiaceous fungi, and yeasts such as *Cryptococcus* spp. and some *Candida* species.⁷
- Some rare non-dermatophyte moulds (*N.dimidiatum*, *N.hyalinum*, *Hortaea werneckii*) are capable of causing dermatophyte-like lesions but are inhibited by cycloheximide. If the clinician mentions the possibility of infection with those moulds, the sample should be plated on a cycloheximide-free medium.⁶
- Chloramphenicol may inhibit some pathogenic fungi.¹⁰
- A single medium is only rarely useful to recover all pathogens contained in a specimen, therefore it is necessary to select media both with and without inhibitory agents for the primary inoculation of the specimen.
- Even if the microbial colonies on the plates are differentiated on the basis of their morphological and chromatic characteristics, it is recommended that biochemical, immunological, molecular, or mass spectrometry testing be performed on isolates, from pure culture, for complete identification. If relevant, perform antimicrobial susceptibility testing.
- This culture medium is intended as an aid in the diagnosis of infectious diseases; the interpretation of the results must be made considering the patient's clinical history, the origin of the sample and the results of other diagnostic tests.

13 - PRECAUTIONS AND WARNINGS

- This product is a qualitative *in vitro* diagnostic, for professional use only; it is to be used by adequately trained and qualified laboratory personnel, observing approved biohazard precautions and aseptic techniques.
- This product is not classified as dangerous according to current European legislation.
- This culture medium contains raw materials of animal origin. The *ante* and *post mortem* controls of the animals and those during the production and distribution cycle of the raw materials, cannot completely guarantee that the product doesn't contain any transmissible pathogen. Therefore, it is recommended that the ready-to-use plates be treated as potentially infectious, and handled observing the usual specific precautions: do not ingest, inhale, or allow to come into contact with skin, eyes, mucous membranes. Download the TSE Statement from the website www.biolifeitaliana.it, describing the measures implemented by Biolife Italiana for the risk reduction linked to infectious animal diseases.
- All laboratory specimens should be considered infectious.
- The laboratory area must be controlled to avoid contaminants such as culture medium or microbial agents.
- Each plate of this culture medium is for single use only.
- Ready-to-use plates are not to be considered a "sterile product" as they are not subject to terminal sterilization, but a product with controlled bio contamination, within the limits of defined specifications reported on the Quality Control Certificate.
- Sterilize all biohazard waste before disposal. Dispose the unused medium and the sterilized plates inoculated with samples or microbial strains in accordance with current local legislation.
- The Certificates of Analysis and the Safety Data Sheet of the product are available on the website www.biolifeitaliana.it.
- The information provided in this document has been defined to the best of our knowledge and ability and represents a guideline for the proper use of the product but without obligation or liability. In all cases existing local laws, regulations and standard procedures must be observed for the examination of samples collected from human and animal organic districts, for environmental samples and for products intended for human or animal consumption. Our information does not relieve our customers from their responsibility for checking the suitability of our product for the intended purpose.

14 - STORAGE CONDITIONS AND SHELF LIFE

Upon receipt, store plates in their original pack at 2-8°C away from direct light. If properly stored, the plates may be used up to the expiration date. Do not use the plates beyond this date. Plates from opened plastic sachet can be used for 7 days when stored in a clean












area at 2-8°C. Do not use the plates if the plastic sachet is damaged or if the dish is broken. Do not use the plates with signs of deterioration (e.g. microbial contamination, dehydration, shrinking or cracking of the medium, atypical colour, excess of moisture).

15 - REFERENCES

1. Espinel-Ingroff A. History of medical mycology in the United States. *Clin Microbiol Rev* 1966;9:235-272
2. Sabouraud R. Contribution à l'étude de la trichophytie humaine. Etude clinique, microscopique et bactériologique sur la pluralité des trichophytons de l'homme. *Ann Dermatol Syphil* 1892; 3:1061-1087.
3. European Pharmacopoeia, current edition.
4. Whiffen AJ, Bonoxos N, Emerson RL. The production of an antifungal antibiotic by *Streptomyces griseus*. *J Bact* 1946; 52: 610-611.
5. Stanley A, Rosenthal D, Furnari BA. The use of cycloheximide-chloramphenicol medium in routine culture of fungi. *J Invest Dermatol* 1957; 28(5):367-71.
6. Public Health England. Investigation of Dermatological Specimens for Superficial Mycoses. UK Standards for Microbiology Investigations. B 39 Issue 3.1, 2016.
7. Lindsley MD. Reagents, stains and media: mycology. In Carrol KC, Pfaller MA et al. editors. *Manual of clinical microbiology*, 12th ed. Washington, DC: American Society for Microbiology; 2019
8. Berkow EL, McGowan KL. Specimen collection, transport and processing: mycology. In Carrol KC, Pfaller MA et al. editors. *Manual of clinical microbiology*, 12th ed. Washington, DC: American Society for Microbiology; 2019.
9. Australian Society for Microbiology: Guidelines for assuring quality of medical microbiological culture media. 2nd Ed, July 2012
10. MacFaddin JF. Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria. Baltimore: Williams & Wilkins; 1985

TABLE OF APPLICABLE SYMBOLS

REF or REF Catalogue number	LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by
 Temperature limitation	 Contents sufficient for <n> tests	 Consult Instructions for Use	 For single use only	 Fragile, handle with care

REVISION HISTORY

Version	Description of changes	Date
Instructions for Use (IFU) - Revision 1	Updated layout and content in compliance with IVDR 2017/746	2021/01
Revision 2	Removal of obsolete classification	2023/03

Note: minor typographical, grammatical, and formatting changes are not included in the revision history.

