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**15 - BIBLIOGRAFIA**

1. Clinical and Laboratory Standards Institute. Reference method for broth dilution antifungal susceptibility testing of yeasts. Approved standard CLSI document M27-A3. Clinical and Laboratory Standards Institute, Wayne, Pa. 2008.
2. Arendrup MC *et al.* EUCAST definitive document E.7.3.2 Method for the determination of broth dilution minimum inhibitory concentrations of antifungal agents for yeasts. Valid from 22 April, 2020.
3. Guinea J *et al.* How to: EUCAST recommendations on the screening procedure E.Def 10.1 for the detection of azole resistance in *Aspergillus fumigatus* isolates using four-well azole-containing agar plates. Clin Microbiol Infect 2019; 25(6):681-687.
4. Arendrup MC *et al.* How to: perform antifungal susceptibility testing of microconidia-forming dermatophytes following the new reference EUCAST method E.Def 11.0, exemplified by *Trichophyton*. Clin Microbiol Infect 2021;27(1):55-60.
5. Ranque SL, Lachaud LM, Gari-oussaint L, Michel-Nguyen A, Mallié M, Gaudart J, Bertout S. Interlaboratory Reproducibility of Etest Amphotericin B and Caspofungin Yeast Susceptibility Testing and Comparison with the CLSI Method. J Clin Microbiol 2012; 50(7): 2305–2309. .
6. Moore GE, Gerner RE, Franklin HA. Culture of normal human leukocytes. JAMA 1967; 199 (8): 519–524
7. Pfaller MA, Messer SA, Mills K, Bolmström A, Jones RN. Evaluation of Etest method for determining caspofungin (MK-0991) susceptibilities of 726 clinical isolates of *Candida* species. J Clin Microbiol 2001;39 (12):4387-9.
8. Espinel-Ingroff A. . Etest for antifungal susceptibility testing of yeasts. Diagn Microbiol Infect Dis 1994;19(4):217-20
9. Pfaller MA, Messer SA, Houston A, Mills K, Bolmström A, Jones RN. Evaluation of the Etest method for determining voriconazole susceptibilities of 312 clinical isolates of *Candida* species by using three different agar media. J Clin Microbiol . 2000;38(10):3715-7.
10. Maxwell MJ, Messer SA, Hollis RJ, Diekema DJ, Pfaller MA. Evaluation of Etest Method for Determining Voriconazole and Amphotericin B MICs for 162 Clinical Isolates of *Cryptococcus neoformans*. J Clin Microbiol. 2003; 41(1): 97–99.
11. Jong Hee Shin JH *et al.* Detection of amphotericin B resistance in *Candida haemulonii* and closely related species by use of the Etest, Vitek-2 yeast susceptibility system, and CLSI and EUCAST broth microdilution methods. J Clin Microbiol 2012;50(6):1852-5.
12. Johnson EM, Cavling-Arendrup M. Susceptibility Test Methods: Yeasts and Filamentous Fungi. In: Carrol KC, Pfaller MA *et al.* editors. Manual of clinical microbiology, 12th ed. Washington, DC: American Society for Microbiology; 2019.
13. Shields RK, Nguyen MH, Press EG, Updike CL, Clary CJ. Anidulafungin and micafungin MIC breakpoints are superior to that of caspofungin for identifying FKS mutant *Candida glabrata* strains and echinocandin resistance. Antimicrob Agents Chemoter 57:6361-6365
14. Axner-Elings M, Botero-Kleiven S, Jensen RH, Arendrup MC. Echinocandin susceptibility testing of *Candida* isolates collected during 1-year period in Sweden. J Clin Microbiol 2011; 49:2516-2521.

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