

EGG YOLK EMULSION 50%

FORMULA (per litre)

Fresh egg yolk500 ml

Saline solution.....500ml

INTENDED USE

Liquid supplement for the differentiation of *Bacillus cereus* or *C.perfringens* . Egg Yolk Emulsion should be used for the supplementation of the following culture media: Bacillus Cereus Agar Base (50 ml/l), Clostridium perfringens Agar Base (100 ml/l). It is used for the detection of lecithinase enzyme activity. See the technical sheets of the above media. For

DIRECTIONS

Bacillus cereus agar MYP: 50 ml Egg Yolk Emulsion + 450 ml Bacillus Cereus Agar Base MYP (REF 401111) autoclaved and cooled to 45-50°C

Bacillus cereus agar PEMBA: 25 ml Egg Yolk Emulsion + 475 ml Bacillus Cereus Agar Base PEMBA (REF 401112) autoclaved and cooled to 45-50°C

SFP Agar: 50 ml Egg Yolk Emulsion + 450 ml Clostridium Perfringens Agar Base (REF 401307) autoclaved and cooled to 45-50°C

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date. Separation may occur on storage. This occurrence is a normal one and doesn't affect the quality of the product. Mix well before the use

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PACKAGING

| | | |
|-----------------|-------------------------------|---------------|
| 42111601 | Egg Yolk Emulsion 50% | 50 ml |
| 42111605 | Egg Yolk Emulsion 50%, | 100 ml |
| 42111600 | Egg Yolk Emulsion 50%, | 200 ml |