Biolife Italiana Srl

IODINE SOLUTION

BOTTLE CONTENTS

lodine	10 g
Potassium Iodide	12.5 g
Distilled Water	50 ml

DIRECTIONS

lodine Solution is ready to use

Intended use

lodine Solution is used to supplement Tetrathionate Broth Base (20 ml of solution/litre of base medium) and Mueller Kauffmann Tetrathionate Broth base (19 ml of solution/litre of base medium). For the details of the procedures see the technical sheets of Tetrathionate Broth Base (cat. N° 402125) and Mueller Kauffmann Tetrathionate Broth Base (cat. N° 401743-401745)

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Packaging

421501 lodine Solution 50 ml