

MUELLER HINTON AGAR BLOOD SHEEP

Ready to use plates

TYPICAL FORMULA * (g/l)

Beef Extract	2.0
Acid Digest of Casein	17.5
Starch	1.5
Agar Bios Special	17.0
Defibrinated sheep blood	50.0 ml

*Adjusted or supplemented to meet performance standard

pH 7.3 ± 0.1

DESCRIPTION

Mueller Hinton Agar Blood Sheep, 90mm and 140mm diameter dishes, are recommended by NCCLS for the antimicrobial disc sensitivity test of streptococci, which are difficult to cultivate on Mueller Hinton Agar without blood.

Consult appropriate references for the procedure of performing antimicrobial disc sensitivity test.

STORAGE

Store at 2-8° - When stored as directed the plates remain stable until the expiry date shown on the label. Do not use beyond stated expiry date.

REFERENCES

• NCCLS Performance Standard for Antimicrobial Susceptibility Testing; Twelfth Informational Supplement. NCCLS Document M100 – S12, January 2002

PACKAGING

541743	Mueller Hinton Agar Blood Sheep	20 ready to use 90mm plates
501743P	Mueller Hinton Agar Blood Sheep	5 ready to use 140mm plates