

SALINE SOLUTION

Ready to use tubes

INTENDED USE

Ready to use Saline Solution for general purposes in microbiological techniques

TYPICAL FORMULA

Sodium Chloride	8.5 g
Distilled water	1000 ml

DESCRIPTION

Saline Solution is available as ready to use sterile tubes. It is recommended in microbiology for the preparation of microbial suspensions, for the dilution of test strains and other applications.

PRECAUTIONS

For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilize all biohazard waste before disposal.

STORAGE

Store at 2-8° away from direct light - When stored as directed the tubed media remain stable until the expiry date shown on the label. Do not use beyond stated expiry date. Media should not be used if there are any signs of deterioration (shrinking, cracking, discoloration) or contamination.

PACKAGING

554610 Saline Solution

20 x 9 ml ready to use tubes